
Sample Prep for Real-Time RT-PCR



Eric Lader, Ph.D elader@ambion.com

Marketing Product Manager

Manager, Product Applications Group

Agenda

Sample prep considerations

- Tissue handling: do's and don'ts
- RNA isolation: what's important
- My favorite RNA isolation method
- Genomic DNA contamination

RT considerations you've never considered

- Endogenous priming
- RT inhibition of PCR
- RT reaction times

Isolating RNA: Issues and Concerns

Performance

suitability for downstream applications

Yield

rapid tissue preservation rapid and thorough sample disruption

Quality


intactness: rapid sample preservation and disruption

purity: removal of contaminants and genomic DNA

stability: complete inactivation of RNase

Isolating RNA: Stability?

Variable Stability of Prepared RNA



Steps in RNA Isolation

- 1. Obtain Sample
 - process immediately
 - preserve it
- 2. Disrupt sample and release RNA
 - use of physical force to break cells (e.g. polytron, grinding)
 - homogenization in chaotrope (GuSCN, LiCl) or detergent (SDS)
- 3. Extract and purify RNA
 - acid phenol/chloroform (RNAwiz, TRIzol, etc.)
 - glass filter (RNAqueous, RNeasy, etc.)
 - oligo (dT) chromatography (PolyAPurist, etc.)
- 4. Store RNA (0.1X TE, 0.1mM EDTA, or 1 mM Citrate pH6.4)

Step 1 in RNA Isolation: sample handling

Obtain tissue sample

Fresh tissue: can not pause until lysate obtained

or ...

Snap Freeze with Liquid N_2 :

Not possible in many settings

Cumbersome; must carry **Dewar**

Grinding tissue is laborious and potentially hazardous

Handling and transfer difficult (cross-contamination?)


Processing large fragments is especially problematic

Step 1 in RNA Isolation: sample handling

Obtain tissue sample

Fresh tissue: can not pause until lysate obtained

(well, not exactly....)



bActin

Mouse Liver 25°C, Processed at indicated time points

Step 1 in RNA Isolation: sample handling

Obtain tissue sample

Fresh tissue: can not pause until lysate obtained

or ...

Snap Freeze with Liquid N_2 :

Not possible in many settings


Cumbersome; must carry Dewar

Grinding tissue is laborious and potentially hazardous

Handling and transfer difficult (cross-contamination?)

Processing large fragments is especially problematic

Effects of tissue freeze/thaw on RNA quality



RNA*later*TM

Tissue preservation and RNA Stabilization Solution

- Aqueous, non-toxic tissue preservation
- Standardizes tissue preservation & nucleic acid isolation
- Provides spatial and temporal separation of collection and processing with no penalties in quality or throughput
- Samples in RNA later can be stored:


at 37°C for one day

at 25°C for 1-2 weeks: ship samples at ambient temperature

at 4°C for months

frozen indefinitely: archival storage of samples

Quality of RNA from RNA later-treated Tissue



RNA Isolation: Sample Handling Summary

Process fresh samples as quickly as possible... (some tissues may be more stable than others!)

.. or preserve in RNAlater and treat as fresh.

Still must be processed quickly into RNAlater, but can then be handled safely.

.. or snap-freeze and process frozen.

Must never thaw: however some small samples can be directly homogenized with a polytron.

Step 2 in RNA Isolation: Sample Disruption

Choice of disruption method critical for yield and quality

Dear Dr. Lader:

I am working at Johns Hopkins School of Medicine. Now we want to isolated RNA from mouse brain cortex for microarray. I have some problems for that and need your troubleshooting and some suggestion. I isolated total RNA following TRIZOL protocol, homogenizing tissue with Sonicator. After redissovling the RNA pellet with DEPC water, RNA concentration was round 2ug/ul. O.D. A260/280 was round 2.0. But when running the RNA at formaldehyde agarose gel (Northernax 10x MOPS gel running buffer, and Northernax formaldehyde load dye were from ambion a half year ago), the smear bands or nothing showed on the gel. I don't know what was wrong during isolation. Before I had worked another lab, my RNA quality was very good for RT-PCR, Northern, and RPA. Only different thing was using Polytron for homogenizing.

Your kind and help are greatly appreciated.

Best regards,

Step 3 in RNA Isolation: RNA recovery

Total RNA Isolation Methods

- Rapid
 - one-step phenol-based scalable
 - glass-binding higher potential throughput
- Phenol-free
 - glass-binding
- Difficult tissues, specific challenges
 - multi-step phenol based

The Hybrid RNA Isolation Method (for people who RNA isolation)

Requires: RNAqueous or equivalent Acid Phenol: Chloroform

Protocol


- Disrupt sample in GiTC lysis buffer with Polytron
- Extract sample with 1:1 acid phenol:chloroform, spin and collect aqueous phase
- Continue with filter-based protocol
- High yields
- Low residual protein (260:280 of ~2.0)
- Vacuum manifold adaptable
- Scalable 1mg 1g tissue !!!
- Handles difficult tissues well
- No filter clogging problems

g DNA Contamination

- Filter-based purification methods yield RNA that is typically 1-10% DNA (based on real-time data).
- There is no RNA isolation method that generates RNA completely free of DNA contamination. Therefore you must DNase I treat your RNA samples
- A follow-up problem is how to get rid of the DNase I and divalent cations so they won't be present during subsequent cDNA synthesis

DNase I treatment of RNA increasing DNase I and time

This sample starts out at 4% DNA



	fold
DNase I/time	<u>reduction</u>
RT+	-
0	0
1X/30'	~100
1X/60'	~250
2X/30'	~250
4X/15'	~100


99.6% removal (0.016% DNA)

G3PDH TaqMan assay

1X DNase = 0.02 units/ul

Divalent Cations and Heat Degrade RNA


- DNase I Buffer with 0.1 mM CaCl₂, 2.5 mM MgCl₂
- Heated 90C, 5 minutes
- Run on formaldehyde agarose gel



- Don't heat kill DNase I
- Unless you remove the divalents, don't heat your RNA in the RT reaction

DNA-free TM **DNase** Treatment and Removal Reagents

- Optimized DNA digestion reagents
- Inactivates DNase without heating, phenol extraction or precipitation
- DNase Removal Reagent is added directly after DNase digestion: inactivates DNase and removes divalent cations (Ca++, Mg++)



Priming 1st strand synthesis

General:

Oligo dT, anchored dT (3' bias, library construction, 3' RACE) Random Primers (non biased distribution)

Gene Specific Primers (more sensitivity?)


Specialized:

Allele specific primers

Functionalized primers (e.g. T7 Promoter + dT, aRNA)

Primers with restriction sites on the ends (cloning)


Efficiency of different 1st strand primers on standard RT reactions



- •random decamers
- •oligo dT
- •G3PDH primers
- •no primers!!

42°C MMLV RT, G3PDH Taqman PCR

Minimizing endogenous priming in RT Reactions



- •MMLV 42°C
- •SSII 48°C
- •AMV 52°C


For standard qRT-PCR, random priming at 42°C is fine.

If there is any reason that you need specific priming:

- a specific site (gene specific or allele specific priming)
- the 3' end (for 3' RACE or other 'anchored' applications)
- using a chimeric or bifunctional primer (T7, restriction site)

...a high stringency RT reaction may be critical.

RT inhibition of real-time PCR




- •No RT
- •AMV (2 U)
- •MMLV (20 U)

•SSII (20 U)

Reverse transcriptase (1/10th the amount in a standard RT reaction) added directly into a 25 ul G3PDH PCR

RT inhibition of real-time PCR:


can be relieved by heat-denaturation



MMLV (20 U)


RT inhibition of real-time PCR

Heat inactivation must be prior to addition of Taq



RT inhibition of real-time PCR:

can be relieved by addition of carrier protein




RT inhibition of real-time PCR

Take home message –

- 1 tube RT-PCR must use carrier protein (RNase inhibitor, αCry)
- 2 tube RT-PCR can either heat kill or use carrier protein

cDNA synthesis: yield vs. time 42°C



Reaction	Relative
time (min)	cDNA
0	0
7.5	30
15	33
30	42
60	51
120	68
240	100


1 ug RNA, 200 units MMLV, random primed β Actin TaqMan assay

cDNA synthesis: yield vs. time 25°C



25°C Time course, MMLV, β Actin TaqMan assay

cDNA synthesis: 25°C vs. 42°C





elader@ambion.com