

PCR ARRAY HANDBOOK v3.0

The Complete Technical Guide to PCR ARRAYS

Cancer
Cytokines
Biomarkers
ECM & Adhesion
Oxidative Stress
Signal Transduction
Inflammation
Stem Cells
MicroRNA
Epigenetics
Toxicology





NEW 100 Selected Peer-reviewed Publications

 **SABiosciences™**

FOCUS ON YOUR PATHWAY™

QUICK REFERENCE

	APOPTOSIS AVAILABLE ARRAYS 27	page 12
	BIOMARKERS AVAILABLE ARRAYS 22	page 12
	CANCER AVAILABLE ARRAYS 28	page 12
	CELL CYCLE AVAILABLE ARRAYS 23	page 12
	CYTOKINES AVAILABLE ARRAYS 22	page 13
	ECM/ADHESION AVAILABLE ARRAYS 25	page 13
	NEUROSCIENCE AVAILABLE ARRAYS 20	page 13
	SIGNAL TRANSDUCTION AVAILABLE ARRAYS 27	page 13
	STEM CELL/DEVELOPMENT AVAILABLE ARRAYS 26	page 13
	TOXICOLOGY/DRUG ADME AVAILABLE ARRAYS 19	page 13

ORDERING INFO

TEL: **301.682.9200**

FAX: **301.682.7300**



www.SABiosciences.com
support@SABiosciences.com

SABiosciences Corporation
6951 Executive Way, Frederick, MD 21703 USA

Compatible Instrument PCR Array Plate Format

Applied Biosystems (ABI)

ABI 5700	A
ABI 7000	A
ABI 7300	A
ABI 7500 Standard 96-well Block	A
ABI 7500 FAST 96-well Block	C
ABI 7700	A
ABI 7900HT Standard 96-well Block	A
ABI 7900HT FAST 96-well Block	C
ABI 7900HT 384-well Block	E
ABI StepOnePlus™	C

Bio-Rad

iCycler	A
iQ5	A
MyiQ	A
CFX96	D
CFX384	E
Chromo4™	A
Opticon 2®	D

Stratagene

Mx3005p	A
Mx3000p	A
Mx4000	D

Roche

LightCycler 480 96-well Block	F
LightCycler 480 384-well Block	G

Eppendorf

Mastercycler ep realplex	A
--------------------------	---

TaKaRa

TP-800	A
--------	---

Fluidigm

BioMark	H
---------	---



PCR ARRAY GUARANTEE

PCR ARRAYS have a truly remarkable performance record. SABiosciences is so confident that you will love the sensitivity, reproducibility, & reliability of our PCR ARRAY System that we back up all of our PCR ARRAY products with a 30 day money-back guarantee. If you are dissatisfied with the performance of your PCR ARRAYS within 30 days of purchase, you are eligible for a full refund.*

* Visit for Complete Details: www.SABiosciences.com/guarantee.php

Pricing is only valid for customers in the USA. Prices are subject to change. Researchers outside the USA should contact their local SABiosciences' distributor for technical information and pricing. The international distributor list can be found online:

<http://www.SABiosciences.com/distributor.php>

PCR ARRAY HANDBOOK VERSION 3.0

RT² PROFILER™ PCR ARRAYS	2
The Most Accurate and Sensitive Technology for Pathway-Powered Gene Expression Analysis	
RT² PROFILER PCR ARRAY SYSTEM	4
A Complete and Reliable Solution for Gene Expression Analysis	
RT² PROFILER PCR ARRAY PERFORMANCE	5
Reliable and Reproducible Pathway-focused Gene Expression Analysis	
NANO PREAMP™ cDNA SYNTHESIS KIT FOR SMALL SAMPLES	6
Enabling the Analysis of Total RNA with as little as 1 ng of RNA	
FFPE PREAMP PCR ARRAY SYSTEM	7
Enabling Pathway-Focused Gene Expression Analysis of FFPE Samples with PCR Arrays	
CUSTOM PCR ARRAYS	8
Focus on Your Genes with PCR Arrays Created from Your Gene List	
CHAMPIONCHIP™ PCR SYSTEM	9
Reliable Chromatin IP with Real-time PCR Precision Achieved in Only One Day	
RT² MICRORNA PCR ARRAYS	10
Simultaneous Detection of Genome-Wide or Pathway-Focused miRNA	
RESEARCH AREA-FOCUSED PCR ARRAY PRODUCT LISTING	12
Search Our Complete Catalog of PCR Array Products by Your Area of Research	
METHYL-PROFILER™ PCR ARRAY SYSTEM	14
Simple, Fast and Reliable DNA Methylation Analysis Without Bisulfite Conversion	
RT² SYBR® GREEN qPCR MASTER MIXES	16
High-Performance SYBR Green Quantitative PCR	
PCR ARRAY DATA ANALYSIS GUIDE	17
Complete Solutions for Interpreting PCR Array Data	
SERVICE CORE FOR PCR ARRAY GENE EXPRESSION ANALYSIS	17
Complete Set of Quick and Convenient Gene Expression Profiling Services	
PCR ARRAY PUBLISHED LITERATURE	
Publications Citing RT ² Profiler PCR ARRAYS from SABiosciences	
ANGIOGENESIS - APOPTOSIS - ATHEROSCLEROSIS - BREAST CANCER - CANCER BIOLOGY	18
CELL CYCLE - CHEMOKINES - COMMON CYTOKINES - CUSTOM PCR ARRAYS	19
DNA DAMAGE - DRUG METABOLISM - ECM - ENDOTHELIAL CELL - GROWTH FACTORS	20
HYPOXIA - INFLAMMATORY CYTOKINES & RECEPTORS	20
INFLAMMATORY RESPONSE & AUTOIMMUNITY - INSULIN - INTERFERON α, β RESPONSE	21
IFN - JAK/STAT - LIPOPROTEIN & CHOLESTEROL - MAPK - NEUROSCIENCE - NFκB	21
OXIDATIVE STRESS - OSTEOGENESIS - P53 - SIGNAL TRANSDUCTION - T CELL ENERGY	22
TGFβ BMP - Th1, Th2, Th3 - Th17 AUTOIMMUNITY/INFLAMMATION - TLR	22
TRANSCRIPTION FACTORS - TOXICOLOGY - TUMOR METASTASIS - WNT SIGNALING	23
COMPLETE PCR ARRAY PRODUCT AND PRICING GUIDE	24
Search Our Complete Catalog of PCR Array Products and Accessories	

RT² PROFILER™ PCR ARRAYS

The Most Accurate and Sensitive Technology for Pathway Gene Expression Analysis



What Are PCR ARRAYS?

RT² Profiler™ PCR Arrays are the most reliable and sensitive gene expression profiling technology for analyzing focused panels of genes in signal transduction, biological process, or disease-related pathways using real-time PCR.

How Are PCR ARRAYS Utilized?

RT² Profiler PCR Arrays are exceptionally valuable tools for studying cancer, immunology, stem cells, toxicology, biomarker discovery/validation, and phenotypic analysis of cell and tissue samples (fresh, frozen, and fixed).

Why Use RT² Profiler™ PCR ARRAYS?

- Simplicity:**
 The simplicity of RT² Profiler PCR Arrays makes routine expression profiling practical in any research laboratory with a real-time instrument.
- Performance:**
 RT² Profiler PCR Arrays have the sensitivity, reproducibility, specificity, and reliability to accurately profile multiple genes simultaneously in 96- or 384-well formats.
- Relevance:**
 RT² Profiler PCR Arrays focus on profiling the genes relevant to the pathways or disease states important to your research.

Anatomy of a 96-well RT² Profiler PCR ARRAY

HUMAN INFLAMMATORY CYTOKINES AND RECEPTORS PCR ARRAY






Figure 1: Each Well in a PCR Array Measures the Expression of a Gene Related to a Pathway or Disease State. Each cataloged PCR Array contains a list of the pathway-focused genes as well as five housekeeping (reference) genes on the array. Wells H6 through H12 contain a panel of proprietary controls to monitor genomic DNA contamination (HGDC) as well as the first strand synthesis (RTC) and real-time PCR efficiency (PPC).

How PCR ARRAYS Work

1. Convert Total RNA to cDNA.



2. Add cDNA to RT² SYBR® Green Master Mix. Aliquot Mixture Across PCR Array.



3. Run in Your Real-Time PCR Instrument.



4. Data Analysis.



SYBR® Green Versus TaqMan® Chemistries




Figure 2: Comparable Biological Results.* Gene expression analysis was compared between RT² Profiler PCR Arrays (SYBR Green-based) and the TaqMan platform. Regression analysis of fold differences, with data normalized against POLR2A, demonstrate that both platforms yield similar biological results.




Figure 3: Sensitivity with RT² SYBR Green Versus TaqMan Chemistry.* PCR amplicons detected using the same primer pair with or without TaqMan probes in either SYBR Green or TaqMan chemistry. SYBR green chemistry yields earlier Ct,s for each dilution, demonstrating better sensitivity than TaqMan chemistry.

* BMC Genomics 2008, 9: 378.

Application: Angiogenesis




Figure 4: Relative Fold Change Between Disorganized and Organized Colonies Using the RT² Profiler Angiogenesis PCR Array. RNA isolated from unorganized T4-2 cells treated with a control antibody (IgG) or reverted to an organized colony by blocking EGFR signaling (mAb225) was reverse transcribed and relative gene expression data was obtained using the Human Angiogenesis PCR Arrays. The expression profile of 84 genes relevant to Angiogenesis as well as 5 housekeeping genes was assayed. Fold change calculations were done using SABiosciences' data analysis software which automatically calculates the fold change in gene expression between the treated and control groups.

Application: Immune Response




Figure 5: Common Cytokine PCR Array Identified 23 Up-Regulated and 6 Down-Regulated Genes Following PBMC Stimulation. Triplicate total RNA samples from human peripheral blood mononuclear cells (either untreated or stimulated with 50 ng/mL PMA and 1 μg/mL ionomycin for 6 hours) were characterized with the Human Common Cytokine PCR Array. Twenty-three cytokine genes are up-regulated (> 5-fold, p < 0.0005) including interleukins, colony stimulating factors, and TNF ligands after 6 hours of stimulation. Six interleukin and TNF ligand genes are down-regulated under the same conditions.

Application: Determining Drug Toxicity with PCR ARRAYS




Figure 6: Stress and Toxicity PathwayFinder™ PCR Array Uncovered Distinct Gene Expression Profiles Associated with Liver Toxicity Caused by 3 PPARγ Agonists. RNA from HepG2 cells treated with three different glitazone PPARγ agonists for type 2 diabetes mellitus was characterized, and the results were compared to that of a vehicle (DMSO) control. The drug withdrawn due to idiosyncratic liver toxicity (Rezulin), induces very different changes in the expression of stress-related genes than two safer drugs still on the market (Avandia and Actos).

Application: ECM PCR ARRAYS for Cancer Biomarker Discovery




Figure 7: ECM and Cell Adhesion PCR Arrays Revealed Up- and Down-Regulated Genes in Breast Cancer. Total RNA from a normal human breast and a human breast tumor were characterized in technical triplicates, and the relative expression levels for each gene in the two samples are plotted against each other in the Scatter Plot. Genes encoding the matrix metalloproteinases (MMP3 and MMP9) and their inhibitors (TIMP3) are up-regulated, while genes encoding integrins (ITGB3 and ITGB4) are down-regulated, by at least three-fold (outside the silver field) in breast tumors relative to normal tissue.


Popular Pathway-Focused PCR ARRAYS	
Angiogenesis	Inflammatory Cytokines & Receptors
Apoptosis	NFκB Signaling Pathway
Cancer PathwayFinder™	Oxidative Stress & Antioxidant Defense
Cell Cycle	Signal Transduction PathwayFinder™
Chemokines and Receptors	Stem Cell
Common Cytokines	Stress and Toxicity PathwayFinder™
DNA Damage Signaling Pathway	TGFβ / BMP Signaling Pathway
Endothelial Cell Biology	Th1-Th2-Th3
Epithelial to Mesenchymal Transition	Toll-Like Receptor Signaling Pathway
Extracellular Matrix and Adhesion	Wnt Signaling Pathway
Custom PCR Arrays <i>(Detailed Information on Page 8)</i>	Gene Expression Analysis Services <i>(Detailed Information on Page 17)</i>

Complete the RT ² Profiler PCR ARRAY System	
Plate Format	Pack Size
96-well Plate	2 Arrays
96-well Plate	12 Arrays
96-well Plate	24 Arrays
384-well Plate	4 Arrays
Add up to 4 Genes to any PCR Array	

PCR ARRAY Accessories	Pack Size	Catalog #
RT ² First Strand cDNA Synthesis Kit	12 Samples	C-03
RT ² Nano PreAMP cDNA Synthesis Kit	12 Samples	C-06
RT ² FFPE PreAMP cDNA Synthesis Kit	12 Samples	C-07
RT ² Nano PreAMP cDNA Synthesis Primer Mixes	12 Samples	Various
SYBR Green w/ ROX Master Mix	2 Arrays	PA-012
SYBR Green w/ Fluorescein Master Mix	2 Arrays	PA-011
SYBR Green Only Master Mix	2 Arrays	PA-010
FREE PCR Array Data Analysis Software		


RT² PROFILER PCR ARRAY SYSTEM

A Complete & Reliable Solution for Gene Expression Analysis



How The PCR ARRAY System Works

SABiosciences RT² Profiler PCR Arrays are a complete system for Pathway-Focused Gene Expression Analysis. From Sample Preparation to Data Analysis, the PCR Array system includes four components that GUARANTEE high-quality, reproducible, and reliable gene expression data.



Integral to the performance of the PCR Array system is a proprietary set of control elements that enhance the reliability of your data and serve as a guarantee for performance over time. These elements allow researchers to quickly assess the quality of their data by determining if samples were contaminated with genomic DNA (gDNA), the quality of the reverse transcription reaction, and real-time PCR efficiency. Each component of the RT² Profiler system contributes to these quality control elements by incorporating an interlocked system for comprehensive monitoring of each step of the PCR Array process.

- **RT² Profiler PCR ARRAYS**

Each pathway-focused PCR Array includes 89-wet bench validated qPCR Primers Assays (including 5 housekeeping genes) and a proprietary control panel.

- **RT² SYBR Green qPCR Master Mixes**

A unique formulation of buffers (GE) that co-evolved with the primer design algorithm provides high amplification efficiencies. Available with reference dyes (ROX, Fluorescein or without).

- **RT² First Strand cDNA Synthesis Kit**

An External RNA Control detected by the PCR Array tests the quality of input RNA. It also features a proprietary genomic DNA elimination buffer essential for eliminating residual gDNA, ensuring specific detection of mRNA.




Figure 1: Elimination of Genomic DNA Contamination. RNA from HEK 293T cells, mouse spinal tissue, mouse brain tissue, or rat brain tissue was characterized on SYBR Green PCR Arrays before (blue bars) and after (red bars) treatment with gDNA Elimination Buffer from the RT² First Strand Kit.





Figure 2: Monitoring Inhibition in Reverse Transcription. Human universal RNA was added with magnesium salt to simulate RNA degradation or added with TRIzol[®] reagent to simulate contamination that inhibits enzyme activity. RT² First Strand Kit was used for cDNA synthesis.

- **FREE Data Analysis Software**

The power of the PCR Array to assess the expression of a pathway-focused set of genes over a wide range of detection yields an abundance of data. With our FREE PCR Array Data Analysis tool, go from raw C_t values to fold change results displayed in a variety of formats (Scatter Plots, Volcano Plots, Cluster-gram) in a **MATTER OF MINUTES**.

RT² PROFILER PCR ARRAY PERFORMANCE

Reliable and Reproducible Pathway-focused Gene Expression Analysis



RT² Profiler PCR Arrays are used and trusted by thousands of research scientists for pathway-focused gene expression analysis. Several factors, including the RT² Primer Assay design algorithm, the proprietary control panel, and the strict manufacturing and quality control procedures, ensure the outstanding performance and reliability of our PCR Arrays. Each PCR Array and every qPCR Primer Assay is wet-bench validated to guarantee their performance, with results demonstrating several performance parameters demonstrated here.

Distinct Specificity

The complete PCR Array System, with high quality input RNA, is guaranteed to yield single bands without primer dimers or other secondary products. The proprietary primer design algorithm incorporates more than ten thermodynamic and sequence alignment criteria, and our wet-bench validation provides confidence that every real-time qPCR Assay accurately represents the expression of the queried gene. Over 20,000 gene-specific RT² PCR Primer Assays have been designed and shipped to satisfied customers.




Figure 1: PCR Arrays Amplify A Single Gene-Specific Product in Every Reaction. Universal total RNA was characterized for four chemokine and chemokine receptors using RT² Primer Assays, followed by a dissociation (melt) curve analysis. PCR Arrays specifically detect individual genes despite the expression of related gene family members in the same RNA sample.

High Sensitivity and Wide Dynamic Range

A key benefit of using pathway-focused PCR Arrays for gene expression analysis is that genes that are over expressed can be measured as reliably as those that are under expressed. The complete PCR Array System yields > 85% positive call with 25 ng - 5 µg RNA or >90% with as little as 1 ng PreAMP RNA. The 8-log wide dynamic range provided by real-time PCR is unparalleled when comparing a pathway-focused gene panel of varying expression levels across a variety of samples.




Figure 2: PCR Arrays Detect as Little as 1 ng RNA. Different amounts of universal total RNA were characterized using the Human Inflammatory Cytokines and Receptors PCR Array (PAHS-011) with or without PreAMP. The percentage of detectable genes was calculated for each RNA amount, with 1ng RNA analysis enabled with the new PreAMP technology.




Figure 3: PCR Arrays Detect RNA Across a Wide Dynamic Range. Ten-fold serial dilutions of Human CHRNA5 were characterized with the respective RT² qPCR Primer Assay.

Uniform PCR Amplification Efficiency

One prerequisite for PCR Array technology is that the amount of template product doubles with every cycle. The more the assays deviate from this ideal, the error in the fold change calculation (ΔC_t) increases exponentially. Only with consistently high amplification efficiencies can PCR Arrays yield meaningful comparison of gene expression levels of all genes simultaneously. The unique combination of SABiosciences' proprietary primer design algorithm and rigorous testing of every primer assay by hand guarantees the high performance of every primer assay on the PCR Arrays.

Superb Reproducibility

Regardless of user or instrument used, the complete PCR Array System demonstrates strong correlations across technical replicates, lots, and instruments with average correlation coefficients > 0.99 insuring reliable detection of differences in expression between biological samples.




Figure 4: PCR Arrays Yield Highly Reproducible Results. Four replicate sets of raw threshold data (1-4) obtained by two different scientists (A & B) at two different times on Human Drug Metabolism RT² Profiler PCR Arrays are directly compared. The results demonstrate a high degree of correlation ($R^2 > 0.990$).

RT² Profiler PCR Arrays: A Trusted & Reliable System

PCR Arrays have been used by thousands of researchers who have successfully submitted and published their PCR Array results in very high impact journals, including Science, PNAS, Cancer Research, the Nature & Cell family of journals, and others (See Pages 18-23).

RT² NANO PREAMP™ cDNA SYNTHESIS KIT

Enabling the Analysis of 1 Nanogram of RNA with RT² Profiler™ PCR ARRAYS



What Is the RT² Nano PreAMP™ cDNA Synthesis Kit?

RT² Nano PreAMP cDNA Synthesis Kit and Primer Mixes are a breakthrough technology enabling expression analysis starting from as little as 1 ng of total RNA. It employs a proprietary preamplification process to faithfully increase the amount of targeted cDNA for PCR Array analysis. This technology empowers RT² Profiler PCR Arrays to accurately analyze **nanogram** levels of total RNA.

Samples that can NOW be characterized with real-time PCR Arrays include:

- **Laser Captured Microdissection Samples (LCM)**
- **Fine Needle Aspiration Biopsies (FNAB)**
- **Stem Cell Clusters or Embryoid Bodies**
- **Flow Cytometry / Fluorescent-Activated Cell Sorting (FACS)**

Combined with PCR Arrays, the RT² Nano PreAMP cDNA Synthesis Kit and Primer Mixes extends the PCR Array System to accurately analyze a pathway-focused set of genes with as little as 1 ng of total RNA.

- **RT² Nano PreAMP cDNA Synthesis Kit:** Proprietary kits include optimized reagents for first strand cDNA synthesis and preamplification from only 1 ng of total RNA.
- **RT² Nano PreAMP cDNA Synthesis Primer Mixes:** Ready-to-use primer mixes for amplifying pathway-specific cDNA templates on corresponding RT² Profiler PCR Arrays.

Benefits of RT² Nano PreAMP cDNA Synthesis Technology

- **Robust Performance on Small Samples:** Analyze up to 4 different PCR Arrays starting with as little as 1 ng of Total RNA.
- **Easy Workflow and Designed for Routine Use:** Simple and quick procedures with minimal hands-on time to preamplify target templates in under 2 hrs.
- **Superior Sensitivity:** Maximally enhances the sensitivity of RT² Profiler PCR Arrays to analyze limited amounts of RNA.





Figure 1. RT² Nano Preamplification Significantly Increases Sensitivity of Detection from 1 ng RNA Samples. Different amounts of human universal RNA were converted to cDNA with (red) or without (blue) RT² Nano preamplification. The unamplified and preamplified samples were then analyzed on the Human Inflammatory Cytokines and Receptors PCR Array (PAHS-011), which contains 84 pathway-specific assays, plus controls, including 5 assays for housekeeping genes. Threshold cycle values (C_t) were obtained and any genes with a C_t < 35 were considered to be present. Results indicate that with Nano preamplification, a 33.7% increase in positive call rate is observed in samples with as little as 1 ng RNA.

How RT² Nano PreAMP cDNA Synthesis Kits Work

Two Simple Steps:



1. cDNA First Strand Synthesis

This kit provides enough reagents for synthesizing first strand cDNA from 12 different samples of as little as 1 ng total RNA.

2. Preamplification of cDNA for Pathway Specific Genes

Each first strand cDNA synthesis reaction can be amplified by 4 different sets of PCR Array-specific Primer Mixes, allowing gene expression analysis for one sample on as many as 4 different PCR Arrays.

UNBIASED AMPLIFICATION PROCESS & COMPARISON OF ΔC_t VALUES BETWEEN PREAMPLIFIED AND UNAMPLIFIED cDNA




Figure 2. Unbiased Amplification Process - Highly Comparable ΔC_t Values Between Preamplified and Unamplified cDNA from Human Liver Tumor RNA. First strand cDNA was synthesized from 5 ng of human liver tumor RNA. One-quarter of each RT product was used for preamplification with the RT² Nano PreAMP cDNA Master Mix Kit plus the Human Cancer PathwayFinder™ Nano PreAMP Primer Mix. Unamplified cDNA synthesized from 500 ng of the same liver tumor RNA sample was used as the control. PreAMP amplified and unamplified cDNA samples were then analyzed on the Human Cancer PathwayFinder™ PCR Array, and the threshold cycle values (C_t) were obtained. The ΔC_t value for each gene was calculated by subtracting the average C_t value of the five reference genes (B2M, HPRT1, RPL13A, GAPDH, and ACTB) on the PCR Array from the C_t value of each gene of interest. The concordance of the ΔC_t values between preamplified and unamplified samples was evaluated by regression analysis. Data points with C_t ≥ 35 were considered to be absent genes and were excluded from the analysis. The dashed line represents the ideal slope of 1.0. The solid line shows a linear regression fit with the R² and slope indicated. The high correlations between preamplified and unamplified cDNA were also obtained from universal RNA samples (*data not shown*).

RT² Nano PreAMP cDNA Synthesis Products

Product	Catalog #
RT ² Nano PreAMP cDNA Synthesis Kit	C-06
RT ² Nano PreAMP cDNA Synthesis Primer Mixes*	
Human Cancer PathwayFinder™	PBH-2033
Human Apoptosis	PBH-7012
Human Angiogenesis	PBH-4024
Human Mesenchymal Stem Cell	PBH-4082
Human Embryonic Stem Cell	PBH-1081
Human Extracellular Matrix and Adhesion Molecules	PBH-0013
Human Inflammatory Cytokines and Receptors	PBH-4011
Mouse Inflammatory Cytokines and Receptors	PBM-9011
Human Toll-like Receptor Signaling Pathway	PBH-5018
Mouse Toll-like Receptor Signaling Pathway	PBM-0018
RT² Nano PreAMP Primer Mixes for All PCR Arrays	PBX-####

* PreAMP Primer Mixes for Custom PCR Arrays available

RT² FFPE PreAMP™ PCR ARRAY SYSTEM

Enabling the Analysis of FFPE Samples with PCR ARRAYS



Gene Expression Analysis from FFPE Samples

An innovative solution enabling the accurate qRT-PCR analysis of Formalin-fixed Paraffin-embedded (FFPE) samples. The RT² FFPE PreAMP technology utilizes multiplex tandem PCR to preamplify gene-specific cDNA with minimal bias. This kit is intended for preamplification of first-strand cDNA from fragmented total RNA from FFPE samples for gene expression analysis with RT² Profiler PCR Arrays.

The combination of a simplified **Xylene-Free** RNA extraction and a high-fidelity amplification process maximizes recovery of RNA and microRNA (miRNA). RT² Profiler PCR Arrays facilitate easy and reliable expression analysis of genes associated with a biological pathway or a diseased state from FFPE samples.

Benefits of RT² FFPE PreAMP PCR ARRAY System

- **Quick and Efficient:** High quality and high-yield total RNA and miRNA isolation from FFPE samples in 70 minutes
- **Superior Sensitivity:** PreAMP protocol significantly enhances qRT-PCR detection sensitivity for FFPE samples
- **Easy Workflow:** Simple **Xylene-Free** procedure and robust performance

RT² FFPE PREAMPLIFICATION FAITHFULLY REPRESENTS BIOLOGICAL CHANGES





Figure 1. Highly Comparable Gene Expression Fold Change Results between FFPE Preamplified and Unamplified Samples. RNA extracted from FFPE spleen and intestine samples were extracted using the RT² FFPE RNA Extraction Kit and converted to cDNA with and without preamplification. All four cDNAs were analyzed on the Human Cancer PathwayFinder™ PCR Array. The $\Delta\Delta C_t$ comparison and genes with raw C_t values lower than 33 in both unamplified spleen and intestine samples are presented.

RT² FFPE PreAMP Performance

- **Increased positive call rate from FFPE samples**
- **Increased detection of genes previously classified as "Absent"**
- **Unbiased amplification of preamplified genes**
- **Faithful conservation of biological changes**

How the RT² FFPE PreAMP PCR System Works



RT² FFPE PreAMP INCREASES DETECTION OF GENES PREVIOUSLY CLASSIFIED AS "ABSENT"





Figure 2. Genes extracted from FFPE samples previously classified as "Absent" are now detectable after RT² FFPE Preamplification. RNA was extracted from FFPE spleen sample (human) with the RT² FFPE RNA Extraction Kit and reverse transcribed to cDNA using RT² FFPE preamplification (red bars) and without PreAMP (blue bars). Results of the Human Cancer PathwayFinder PCR Array showed 55% of unamplified genes were virtually undetectable with no genes in the 10-20 C_t range. Preamplified genes with C_t values > 30, shift into the reliably quantitative range (C_t = 10-30).

RT² FFPE PreAMP & RNA Extraction Products

Product	Catalog #
RT ² FFPE RNA Extraction Kit	PA-023
RT ² FFPE PreAMP cDNA Synthesis Kit	C-07
RT ² FFPE PreAMP Primer Mixes for All PCR ARRAYS	PFX-####
Custom RT ² FFPE PreAMP Primer Mixes for PCR ARRAYS	Inquire
Pathway-Focused RT ² Profiler PCR ARRAYS	Inquire
Instrument-Specific SYBR® Green qPCR Master Mixes	Inquire

CUSTOM PCR ARRAYS & ASSAY DESIGN

Focus on Your Genes with PCR Arrays & Assays Created from Your Gene List



What Are Custom PCR ARRAYS?

Custom RT² Profiler™ PCR Arrays are a high-throughput approach for profiling the expression of your genes of interest. Choose from any gene in the human, mouse, rat, Rhesus Macaque or Drosophila genomes (up to 384 different genes). Whether your interests are in biomarker discovery, microarray followup, drug development, disease characterization, or signal transduction mechanisms, Custom RT² Profiler PCR Arrays enable focused expression analysis on your genes of interest.

Why Custom PCR ARRAYS from SABiosciences?

- **Performance:** Each assay within a Custom RT² Profiler PCR Array is designed and wet bench-validated using a set of rigorous parameters to insure the genes in your sample across a wide dynamic range are reproducibly recognized and quantified.
- **Flexibility:** Custom RT² Profiler PCR Arrays are available in a number of easy-to-use formats for quick sample loading and data analysis.
- **Turnaround Time:** Submit your 96- or 384-gene list and receive your Custom PCR Arrays in 2 weeks.

Gene Layouts

96-well Custom PCR Arrays

384-well Custom PCR Arrays

Customer Data: Validate Microarray Expression Analyses




Figure 1: Gene Expression Changes in Hypothalamus of MECP2 Mouse Models. Validation of expression changes for 66 genes by qPCR analysis. Gene expression levels from microarray analyses were validated in four MECP2-Tg males and four MeCP2-null males. Data is plotted as relative up-regulation (red) or down-regulation (blue) over Wild-type ($P < 0.05$, t test). Each column represents a single gene, and represents data from four samples for each genotype.

Chahrouh, M., et al. *Science* 2008; 320: 1224-1229.

96-well Custom PCR ARRAYS	
Format	# of Plates (minimum)
12 Genes, 8 Samples / Plate	12, 24
24 Genes, 4 Samples / Plate	24
32 Genes, 3 Samples / Plate	
48 Genes, 2 Samples / Plate	
96 Genes, 1 Samples / Plate	
All Formats	>24

384-well Custom PCR ARRAYS	
Format	# of Plates (minimum)
16 Genes, 24 Samples / Plate	6, 12
32 Genes, 12 Samples / Plate	6, 12
48 Genes, 8 Samples / Plate	
96 Genes, 4 Samples / Plate	
Up to 384 Genes, 1 Samples / Plate	24
All Formats	>24

PCR ARRAY Accessories	Pack Size	Catalog #
RT ² First Strand cDNA Synthesis Kit	12 Samples	CA-03
RT ² qPCR-Grade RNA Isolation Kit	12 Samples	PA-001
RT ² Nano PreAMP cDNA Synthesis Kit	12 Samples	C-06
RT ² FFPE RNA Extraction Kit	12 Samples	PA-023
RT ² FFPE PreAMP cDNA Synthesis Kit	12 Samples	C-07
Custom RT ² Nano or FFPE PreAMP Primer Mixes		Inquire
SYBR [®] Green Master Mixes (see page 16)		
FREE PCR Array Data Analysis Software		

Data shown is used with permission of *Science* and respective authors.

CHAMPIONCHIP™ qPCR SYSTEM

Reliable Chromatin IP with Real-Time PCR Precision Achieved in Only ONE DAY!



The ChampionCHIP™ System provides the first complete platform for the analysis of *in vivo* protein-DNA interactions using Chromatin Immunoprecipitation (ChIP) and real-time PCR (qPCR) detection. With validated high-quality ChIP-grade antibodies and PCR primers for any genes' promoter region, a simple and robust one-day preparation assay quickly delivers reliable and biologically relevant results. ChIP qPCR is a powerful and versatile method for the analysis of chromatin DNA bound by transcription factors, co-regulators, modified histones, chromatin remodeling proteins, or other nuclear factors from live cells. However, the tedious process and variable results have limited many researchers' ability to adopt this technique to study dynamic protein-DNA interactions in native chromatin environments. The ChampionCHIP™ System, yields the most reliable ChIP assay results with qPCR precision in just a single day.

How the ChampionCHIP™ System Works

The ChampionCHIP One-Day Kit simplifies the usual two- to five-day ChIP protocol down to a manageable six to eight hours. Its crosslink reversal step is much faster and less tedious than conventional methods, and its DNA purification step yields a larger quantity and higher quality ChIP DNA than other one-day kits.




Figure 1: The Entire ChampionCHIP System Protocol Can Be Completed in a Single Day. The ChampionCHIP System includes a simplified high-performance One-Day ChIP Kit, ChIP-Grade Antibodies, real-time PCR primers, and a FREE ChIP PCR Array Data Analysis Suite.

High Specificity




Figure 2. The ChampionCHIP System Readily and Correctly Identifies Different Euchromatin and Heterochromatin Loci. ChampionCHIP antibodies against modified histones (H3K4me2, H3K27me3, H3K9me3) or control IgG were used for precipitating chromatin from HeLa cells. Each ChIP DNA fraction was analyzed by real-time PCR using primers specific for the ALDOA, MYO-D, and SAT2 loci to calculate percentages of co-precipitating DNA relative to input.




Figure 3: Differential Histone Modification Tiling Across the Sequence of Any Gene. ChampionCHIP antibodies for modified histones (H3Ac, H3K4me2, H3K27me3) or control IgG were used for precipitating chromatin from one million HeLa cells. Each ChIP DNA fraction was analyzed with a ChampionCHIP Tiling Array representing 30 one-kb tile intervals across the genomic sequence of the CDKN1A gene. The results obtained from three independent experiments are consistent with active transcription of the CDKN1A gene.




Figure 4: Treatment with 5-Fluorouracil Increases CDKN1A Gene Expression and p53 Binding in Cell Lines Expressing Wild-Type But Not Mutant p53. Triplicate samples from A549, HepG2, and PC3 cells were treated with 5-FU (300 μM, 6 h), and either subjected to ChIP with an anti-p53 antibody followed by qPCR analysis of the CDKN1A-2kb p53 binding site, or harvested for RNA to analyze CDKN1A expression by real-time RT-PCR. The results of both assays are expressed as the fold-increase upon 5-FU treatment.

ChampionCHIP qPCR System	
Product	Catalog #
ChampionCHIP PCR Arrays	
Human Stem Cell Transcription Genes	GH-501
Human Oncogene & Tumor Suppressor Genes	GH-502
Human T Helper Cell Differentiation	GH-503
ChampionCHIP One-Day Kit	
ChampionCHIP Antibody Kits (Antibody + qPCR Controls)	See Website
Human RNA Polymerase II, p53 and Histones (H3Ac, H4Ac, H3Kme2, H3K9me3, H3K27me3)	
ChampionCHIP qPCR Primers	Inquire
Any Promoter Region in the Human, Mouse or Rat Genome	
Custom ChampionCHIP qPCR Array	Inquire
Transcription Factor Binding Sites or Promoter Tiling	
SYBR Green w/ ROX Master Mix	PA-012
SYBR Green w/ Fluorescein Master Mix	PA-011
SYBR Green Only Master Mix	PA-010
FREE ChIP PCR Array Data Analysis Software	

RT² MICRORNA PCR ARRAYS

Simultaneous Detection of Genome-Wide or Pathway-Focused miRNA



What are miRNAs?

MicroRNAs are endogenous single-stranded RNA molecules 19-25 nucleotides in length, synthesized in a regulated manner from larger RNA molecules. Many miRNA sequences have been found in a variety of species [for example, ~700 miRNAs in humans, and ~500 in mice].

Why Use SABiosciences miRNA PCR ARRAYS?

Detecting every miRNA across the entire genome in a specific and sensitive way is a very technologically challenging task. Many miRNA family members and otherwise distinct miRNA species have very similar sequences. Moreover, other RNA species such as snRNA, tRNA, mRNA, and rRNA can cause non-specific amplification, making the specific analysis of mature miRNA even more problematic. SABiosciences' proprietary miRNA detection technologies enable uniformly high PCR amplification efficiencies, allowing simultaneous detection of miRNA under uniform cycling conditions.

The RT² miRNA PCR Array accurately analyzes the expression of up to 96 or 384 microRNA sequences simultaneously on ANY real-time PCR instrument. SABiosciences' patent-pending miRNA technology integrates a universal-tailing and reverse transcription reaction specific for miRNA with accurate expression level measurement of distinct miRNA sequences that may differ by a single nucleotide base. RT² miRNA PCR Arrays are the most specific and sensitive technology for analyzing genome-wide miRNA expression.

- **Sensitivity:** As little as 0.5 µg total RNA needed
- **Multi-Sequence Flexibility:** Analyze up to 384 sequences simultaneously
- **Simplicity:** As easy as a real-time PCR Array experiment

Why Study miRNA?

MicroRNA represents a new layer of regulation in endogenous gene transcription and translation. Since there are ~700 miRNAs in humans, with each miRNA potentially having hundreds of targets, the majority of genes may be subject to regulation by one or more miRNAs. miRNAs are already being considered as cancer biomarkers, and their importance is being realized in a variety of other research areas, such as differentiation, neurobiology and immunology. There are three major ways to start studying miRNA:

1. Expression Analysis:

The best technology for determining the expression of miRNA is the SABiosciences miRNA PCR Array System.

2. Bioinformatic Prediction:

Identification of miRNAs that potentially regulate your genes of interest is possible with our powerful yet simple bioinformatic algorithm at:


<http://www.SABiosciences.com/miRNAsearch.php>

3. Functional Studies:

The function of individual miRNAs can be identified via miRNA over expression or suppression of miRNA function. [Over 500 miRNAs in their endogenous genomic context are available as expression constructs].

How miRNA PCR ARRAYS Work

As Easy as a Real-time PCR Experiment



Specificity




Figure 1: miRNA qPCR Assays Distinguish Single Nucleotide Mismatches. RT² miRNA PCR Assays and a competitor's assays for miR-99a and miR-100 were used to detect both corresponding synthetic templates, whose sequences differ by only one nucleotide (B). Relative detection of the off-target template is calculated as a percentage of the correct template detection (A). RT² miRNA PCR Assays' proprietary primer design specifically discriminates closely related sequences better than competing assays.

Sensitivity and Dynamic Range




Figure 2: 1000-Fold Higher Sensitivity Than Our Competitors. Samples containing serially diluted synthetic miR-658 template spiked into a constant amount of small RNA that lacks miR-658, were analyzed with miR-658-specific RT² miRNA qPCR Assays and other commercial assays. The resulting threshold cycle values are plotted versus the amount of synthetic template used. RT² miRNA PCR Assays' proprietary reaction formulation detects 1000-fold lower concentrations of miR-658 with a wider linear dynamic range than competing assays.

Accuracy




Figure 3: miRNA PCR Assays Yield the Most Accurate Results. Threshold cycle values from miR-16 and miR-21 specific RT² miRNA qPCR Assays are plotted versus the amount of HEK293 small RNA in a serial dilution series. Serial dilutions of pooled synthetic cDNA templates were similarly used to calculate the amplification efficiencies of 468 assays (average efficiency of 95.37% ± 6.19%). Consistently high amplification efficiencies and sensitivity enable miRNA PCR Arrays to accurately analyze multiple sequences simultaneously using the ΔΔC_t method.

Reproducibility




Figure 4: miRNA PCR Arrays Yield Highly Reproducible Results. Duplicate samples of human brain small RNA were characterized with the Human Genome RT² miRNA PCR Array, and the raw threshold cycle values from each array were plotted against each other. The linear fit of the data with a strong correlation factor (RT² > 0.99) indicates that miRNA PCR Array results can be reliably compared between plates, runs, biological replicates, and samples.

Application: Cancer Research - Colon Cancer Biomarkers




Figure 5: RT² Human Cancer miRNA PCR Array Identifies Potential Colon Cancer Biomarkers. Small RNA isolated from human colon tumor and matched adjacent normal tissue (Biochain) were characterized with PCR Arrays containing assays specific for 88 cancer-related human miRNA sequences. Fold-differences are calculated from raw C_t values normalized to a panel of housekeeping small nuclear RNA. Several miRNA biomarkers are up-regulated in colon cancer.

Application: Stem Cell Research - Osteogenic miRNA Markers





Figure 6: Cell Differentiation and Development miRNA PCR Arrays Identify Sequences Potentially Regulating Osteogenesis. Human adipose tissue-derived mesenchymal stem cells (hMSC) grown in normal or osteogenic differentiation medium for 16 days were stained with Alizarin Red (A) and isolated for relative miRNA expression profiling. The results are shown in the scatter plot (B). Many miRNA sequences exhibit unique time-dependent changes in expression compared to undifferentiated controls (C).

RT ² miRNA PCR ARRAYS		
Product	Human	Mouse
Genome, 384-well	MAH-3100	MAM-3100
Cancer	MAH-102	MAM-102
Cell Differentiation / Development	MAH-103	MAM-103
miFinder™	MAH-001	MAM-001
Custom RT ² miRNA PCR Array	Inquire	Inquire


miRNA PCR ARRAY Accessories	Pack Size	Catalog #
RT ² qPCR-Grade miRNA Isolation Kit	12 Samples	MA-01
RT ² miRNA First Strand Kit	12 Samples	MA-03
SYBR® Green w/ ROX® Master Mix	2 Arrays	PA-012
SYBR® Green w/ Fluorescein Master Mix	2 Arrays	PA-011
SYBR® Green Only Master Mix	2 Arrays	PA-010

FREE miRNA PCR Array Data Analysis Software


APOPTOSIS




BIOMARKERS



CANCER



CELL CYCLE



RT² PROFILER™ PCR ARRAYS

- Technology Overview (page 2)
- Over 100 Pathways Available (page 24)
- Start With as Little as 1 ng of RNA (page 6)
- Analyze FFPE Samples with PCR Arrays (page 7)
- Over 600 Peer-reviewed Publications
- Custom PCR Arrays Available for Human, Mouse, Rat, Rhesus Macaque and Drosophila
- To Learn More, Please Visit the PCR Product Web Page:



www.SABiosciences.com/RTPCR.php

Angiogenesis
 Apoptosis
 Autophagy
 Cancer PathwayFinder™
 Cell Cycle
 DNA Damage Signaling Pathway
 Drug Metabolism: Phase II Enzymes
 Endothelial Cell Biology
 Heat Shock Proteins
 MAP Kinase Signaling Pathway
 NFκB Signaling Pathway
 Oxidative Stress and Antioxidant Defense
 p53 Signaling Pathway
 PI3K-AKT Signaling Pathway
 Signal Transduction PathwayFinder
 Stress and Toxicity PathwayFinder
 TNF Ligand and Receptor
 Unfolded Protein Response*

Breast Cancer and Estrogen Receptor Signaling
 Cell Surface Markers
 Dendritic and Antigen Presenting Cell
 Epigenetic Chromatin Modification Enzymes*
 Epigenetic Chromatin Remodeling Factors*
 Epithelial to Mesenchymal Transition (EMT)*
 Extracellular Matrix and Adhesion Molecules
 Hematopoietic Stem Cells and Hematopoiesis
 Homeobox (HOX) Genes
 Mesenchymal Stem Cell
 Stem Cell
 T-cell and B-cell Activation
 Th1-Th2-Th3

Angiogenesis
 Angiogenic Growth Factors & Angiogenesis Inhibitors
 Apoptosis
 Autophagy
 Breast Cancer and Estrogen Receptor Signaling
 Cancer Drug Resistance and Metabolism
 Cancer PathwayFinder
 Cell Cycle
 Chemokines & Receptors
 DNA Damage Signaling Pathway
 Epithelial to Mesenchymal Transition (EMT)*
 EGF / PDGF Signaling Pathway
 MAP Kinase Signaling Pathway
 Oxidative Stress and Antioxidant Defense
 p53 Signaling Pathway
 Protein Phosphatases*
 PI3K-AKT Signaling Pathway
 TGFβ BMP Signaling Pathway
 Tumor Metastasis
 Wnt Signaling Pathway

Angiogenesis
 Apoptosis
 Autophagy
 Cancer PathwayFinder
 Cell Cycle
 DNA Damage Signaling Pathway
 Epithelial to Mesenchymal Transition (EMT)*
 MAP Kinase Signaling Pathway
 Neurogenesis and Neural Stem Cell
 NFκB Signaling Pathway
 p53 Signaling Pathway
 PI3K-AKT Signaling Pathway
 Signal Transduction PathwayFinder
 Transcription Factors

RT² MicroRNA PCR ARRAYS

- Regulation of Gene Transcription and Translation (page 10)



Cancer
 Cell Differentiation and Development
 miFinder™
 Whole Genome

Cancer
 Cell Differentiation and Development
 miFinder™
 Whole Genome

Cancer
 miFinder™
 Whole Genome

Cancer
 Cell Differentiation and Development
 miFinder™
 Whole Genome

Methyl-Profiler™ PCR ARRAYS

- Accurate Detection of DNA Methylation at CpG Islands Without Bisulfite (page 12)



Breast Cancer
 Gastric Cancer
 Liver Cancer
 Lung Cancer
 Prostate Cancer


Breast Cancer
 Gastric Cancer
 Liver Cancer
 Lung Cancer
 Prostate Cancer

Breast Cancer
 Gastric Cancer
 Liver Cancer
 Lung Cancer
 Prostate Cancer


Breast Cancer
 Gastric Cancer
 Liver Cancer
 Lung Cancer
 Prostate Cancer

* NEW PCR Array

CYTOKINES / INFLAMMATION




ECM / ADHESION




NEUROSCIENCE



SIGNAL TRANSDUCTION



STEM CELL / DEVELOPMENT



TOXICOLOGY / DRUG ADME



Chemokines & Receptors
Common Cytokine
Inflammatory Cytokines and Receptors
Inflammatory Response and Autoimmunity
Interferon and Receptor
Interferon α , β Response
JAK / STAT Signaling Pathway
NF κ B Signaling Pathway
T Cell Anergy & Immune Tolerance
T-cell and B-cell Activation
TGF β BMP Signaling Pathway
Th17 for Autoimmunity and Inflammation
Th1-Th2-Th3
Toll-Like Receptor Signaling Pathway
TNF Ligand and Receptor
Unfolded Protein Response*

Angiogenesis
Angiogenic Growth Factors & Angiogenesis Inhibitors
Atherosclerosis
Chemokines and Receptors
Common Cytokine
Embryonic Stem Cells
Endothelial Cell Biology
Extracellular Matrix and Adhesion Molecules
Growth Factors
Inflammatory Cytokines and Receptors
MAP Kinase Signaling Pathway
Mesenchymal Stem Cell
NF κ B Signaling Pathway
Osteogenesis
TGF β BMP Signaling Pathway
Tumor Metastasis
TNF Ligand and Receptor

Alzheimer's Disease
Apoptosis
cAMP / Ca²⁺ Signaling PathwayFinder
Drug Transporters
Embryonic Stem Cells
GPCR Signaling PathwayFinder
Heat Shock Proteins
Hedgehog Signaling Pathway
Hypoxia Signaling Pathway
Mesenchymal Stem Cell
Neurogenesis and Neural Stem Cell
Neuroscience Ion Channels and Transporters
Neurotransmitter Receptors and Regulators
Neurotrophin and Receptors
Nitric Oxide Signaling Pathway
Notch Signaling Pathway
Oxidative Stress and Antioxidant Defense
Stem Cell

cAMP / Ca²⁺ Signaling PathwayFinder™
EGF / PDGF Signaling Pathway
GPCR Signaling PathwayFinder™
Hedgehog Signaling Pathway
Insulin Signaling Pathway
JAK / STAT Signaling Pathway
Lipoprotein Signaling and Cholesterol Metabolism
MAP Kinase Signaling Pathway
NF κ B Signaling Pathway
Notch Signaling Pathway
Nuclear Receptors and Coregulators
PI3K-AKT Signaling Pathway
Signal Transduction PathwayFinder
TGF β BMP Signaling Pathway
Toll-Like Receptor Signaling Pathway
Transcription Factors
Ubiquitination Pathway
Wnt Signaling Pathway

Cell Cycle
Dendritic and Antigen Presenting Cell
Embryonic Stem Cells
Extracellular Matrix and Adhesion Molecules
GPCR Signaling PathwayFinder
Hedgehog Signaling Pathway
Hematopoietic Stem Cells and Hematopoiesis
Homeobox (*HOX*) Genes
Lipoprotein Signaling and Cholesterol Metabolism
Mesenchymal Stem Cell
Neurogenesis and Neural Stem Cell
Neurotrophin & Receptors
Notch Signaling Pathway
Osteogenesis
Stem Cell
T-cell and B-cell Activation
TGF β BMP Signaling Pathway
Toll-Like Receptor Signaling Pathway
Transcription Factors
Wnt Signaling Pathway

Apoptosis
Autophagy
Cancer Drug Resistance and Metabolism
Cancer PathwayFinder™
Cell Cycle
DNA Damage Signaling Pathway
Drug Metabolism
Drug Metabolism: Phase I Enzymes
Drug Metabolism: Phase II Enzymes
Drug Transporters
GPCR Signaling PathwayFinder
Lipoprotein Signaling & Cholesterol Metabolism
Mitochondria
Oxidative Stress and Antioxidant Defense
p53 Signaling Pathway
PI3K-AKT Signaling Pathway
Stress and Toxicity PathwayFinder

miFinder™
Whole Genome

Cancer
miFinder™
Whole Genome

miFinder™
Whole Genome

Cancer
Cell Differentiation and Development
miFinder™
Whole Genome

Cancer
Cell Differentiation and Development
miFinder™
Whole Genome

miFinder™
Whole Genome

Inflammatory Response
T Cell Activation
Cytokine Production


Breast Cancer
Gastric Cancer
Liver Cancer
Lung Cancer
Prostate Cancer

Breast Cancer
Gastric Cancer
Liver Cancer
Lung Cancer
Prostate Cancer

Stem Cell Transcription Factors


METHYL-PROFILER™ PCR ARRAY SYSTEM

Simple, Fast and Reliable DNA Methylation Analysis **WITHOUT BISULFITE!**



The Methyl-Profiler™ DNA Methylation PCR Array System is an innovative technology enabling fast and accurate detection of DNA methylation status at CpG islands. This technology replaces the tedious and inefficient bisulfite-based methods with simple and simultaneous selective restriction digests of either methylated or unmethylated DNA, and takes advantage of the quantitative power of real-time PCR. The PCR Array format reveals the DNA methylation status of gene panels related to diseases or pathways. The individual primer pairs allow analysis of the DNA methylation status of any human or mouse gene. The reliability and simplicity of the procedure makes this technology ideal for profiling DNA methylation and biomarkers of stem cell growth and differentiation, cancer, and other human diseases.

How DNA Methylation PCR ARRAYS Work



The human genome contains many long hypermethylated stretches of CpG dinucleotide-rich sequences. In this sea of CpG methylation, unmethylated CpG-rich sequences, known as “CpG islands”, are found in the promoters of most transcriptionally active genes. These normal patterns of DNA methylation are perturbed in cancer cells, where specific tumor suppressor genes (TSG) become hypermethylated, causing their expression to be silenced. Since every tumor type has a unique “methylation profile”, or panel of hypermethylated genes, the analysis of TSG hypermethylation has become very important for basic cancer research, clinical diagnostics, and therapeutic applications.

The Methyl-Profiler PCR Array System fulfills the need to rapidly and simultaneously determine the methylation status of more genes in more samples in a higher-throughput fashion. The current time- and labor-intensive methodologies require bisulfite conversion of unmethylated cytosines to uracil followed by either sequence analysis or PCR using primers sensitive to the resulting base conversion. Bisulfite conversion is not only tedious but is also inefficient and damages DNA. The resulting low yields of DNA make bisulfite-based methods unsuitable for the analysis of small samples.

Why Use the Methyl-Profiler™ PCR ARRAY System?

- **Simple, Fast and Reliable:**
No bisulfite conversion and ready-to-use.
- **Disease- or Pathway-Focused Gene Sets:**
Simultaneously detect DNA methylation of 24 or 96 genes.
- **Genome-wide Coverage:**
Primers to detect methylation of your favorite genes.

The Methyl-Profiler PCR Array System relies on the differential cleavage of target sequences by two different restriction endonucleases whose activities require either the presence or absence of methylated cytosines in their respective recognition sequences. As real-time PCR quantifies the relative amount of intact DNA remaining after each enzyme digestion, the methylation status of individual genes and the methylation profile across a gene panel are reliably and easily calculated. The high yield of DNA from the restriction digests and PCR amplification allow the analysis of smaller, more heterogeneous samples.




Figure 1: Methyl-Profiler PCR Assays Yield Results Consistent with Bisulfite Sequencing. The methylation status of the cadherin 1 gene (CDH1) was determined using bisulfite sequencing and Methyl-Profiler PCR Assays in three breast cancer cell lines known to have very different CDH1 methylation patterns.

To validate the accuracy of the Methyl-Profiler PCR Array System, its results were compared with those generated by bisulfite sequencing, the gold standard for DNA methylation analysis. The methylation status for both the CDH1 (Figure 1) and CDH13 (*data not shown*) genes observed in three different breast cancer cell lines by the two methods match very closely. Real-time PCR characterization of methylation-dependent and methylation-sensitive restriction enzyme digests directly quantifies unmethylated and hypermethylated genomic DNA, respectively. The results indicate that the sensitivity and specificity of the Methyl-Profiler PCR Array System rivals bisulfite sequencing, suggesting that it can readily replace more tedious bisulfite PCR validation methods.

Make your DNA Methylation analysis as quick and painless as possible with the Methyl-Profiler DNA Methylation PCR System and EXCEL-based data analysis.

Download our **FREE EXCEL** Data Analysis Software:

http://www.SABiosciences.com/dna_methylation.php




Figure 2: Methyl-Profiler PCR Assays Detect Hypermethylation in Heterogeneous Samples Containing As Little As Five Percent Tumor DNA. SKBR3 breast cancer cell line and normal blood genomic DNA (encoding hypermethylated and unmethylated HIC1, respectively) were mixed in different ratios. Using Human HIC1 Methyl-Profiler qPCR Primers, the percentage of hypermethylated HIC1 relative to total promoter DNA in each mixture was detectable down to five percent.

Primary tumors are typically very heterogeneous, containing a mixture of both cancerous and noncancerous cells. Therefore, reliable tumor characterization requires detecting smaller amounts of hypermethylated DNA diluted in an unmethylated background. Methyl-Profiler PCR Assays have the sensitivity required to detect hypermethylated DNA from breast cancer cells even when they represent only five percent of the total cell population (Figure 2).

Applications




Figure 3: Methyl-Profiler PCR Arrays Validate Breast Cancer Gene Methylation Status in Breast Cancer Cell Lines. Heat map comparison of the hypermethylation status of 24 genes in the genomic DNA of three breast cancer cell lines and blood genomic DNA as determined by Human Breast Cancer Signature Panel DNA Methylation PCR Arrays.

To demonstrate that Methyl-Profiler PCR Arrays can validate methylation biomarkers, we first scanned published results to design a cataloged PCR Array representing a signature panel of the 24 most frequently methylated genes in human breast tumors. We then analyzed the methylation profile of this gene panel in three different breast cancer cell lines (Figure 3). The results further strengthen the correlation of these biomarkers with breast cancer.




Figure 4: Methyl-Profiler PCR Arrays Discover New Candidate Breast Cancer DNA Methylation Biomarkers. Heat map comparison of the hypermethylation status of a panel of 79 transcription factor genes in six breast cancer cell lines and a normal epithelial cell line as determined with Custom DNA Methylation PCR Arrays.

To demonstrate that Methyl-Profiler PCR Arrays can also discover new biomarkers, we arranged a custom array containing a panel of candidate transcription factor genes, whose methylation status had not been previously associated with breast cancer (Figure 4). We found that breast cancer cell lines also hypermethylate this gene panel, potentially providing a new discovery source for cancer biomarkers.

The Methyl-Profiler DNA Methylation PCR Array System is ideally suited to genomic DNA hypermethylation analysis for both basic research applications and clinical biomarker development. The simple two-step procedure is considerably faster and easier than current bisulfite sequencing and bisulfite PCR methods, and yields closely matching results with equivalent sensitivity.

Methyl-Profiler™ DNA Methylation PCR ARRAYS	
Product*	Catalog #
Human Breast Cancer - Signature Panel	MeAH-011
Human Gastric Cancer - Signature Panel	MeAH-021
Human Liver Cancer - Signature Panel	MeAH-031
Human Lung Cancer - Signature Panel	MeAH-041
Human Prostate Cancer - Signature Panel	MeAH-051
Human Colon Cancer - Signature Panel	MeAH-061
Human Stem Cell Transcription Factors - Signature	MeAH-511
Human Inflammatory Response - Signature Panel	MeAH-521
Human T Cell Activation - Signature Panel	MeAH-531
Human Cytokine Production - Signature Panel	MeAH-541
Custom Methyl-Profiler PCR Arrays	Inquire

* Methyl-Profiler PCR Arrays are available in Signature Panels (24 genes) & Comprehensive Panels (96 genes).


PCR ARRAY Accessories	Pack Size	Catalog #
Methyl-Profiler qPCR Assays	200 Reactions	Inquire
Methyl-Profiler Enzyme Kit	12 Samples	MeA-03
SYBR Green w/ ROX Master Mix	2 Arrays	PA-012
SYBR Green w/ Fluorescein Master Mix	2 Arrays	PA-011
SYBR Green Only Master Mix	2 Arrays	PA-010
FREE Methyl-Profiler PCR Array Data Analysis Software		

RT² SYBR[®] GREEN qPCR MASTER MIXES

From the Experts in High-Performance SYBR[®] Green Quantitative PCR



SYBR Green Detection is a popular approach used in quantifying gene expression analysis with RT-PCR. It relies on the preferential binding of the SYBR green dye to double-stranded DNA, resulting in strong fluorescence emission signals, with the signal intensity proportional to the amount of double-stranded DNA present.



High quality PCR reaction components are essential for achieving superior amplification specificity and efficiency. SABiosciences offers a complete solution for using SYBR Green PCR Arrays with the RT² SYBR Green qPCR Master Mixes. Each mix includes a Hot Start Taq DNA polymerase, which provides tighter control over activity, and other proprietary chemical components that significantly minimize primer dimer formation, thereby enhancing amplification efficiencies for even the most difficult-to-amplify genes. The higher SYBR Green signal from our formulations provides greater sensitivity and ensures clean results without sacrificing specificity or amplification efficiency.

Brighter SYBR Green Signal




Figure 1: RT² SYBR Green qPCR Master Mixes Provide Greater Sensitivity with a Brighter SYBR Green Signal. Four commercial master mixes were used to detect the expression of human ACTB from the same universal reference RNA. The amplification (A) and the dissociation curves (B) for the master mix from SABiosciences demonstrate a sharper amplification curve and a brighter SYBR Green signal than observed with three competing master mixes.

Greater Sensitivity Without Sacrificing Specificity




Figure 2: RT² SYBR Green qPCR Master Mixes Provide Greater Sensitivity Without Sacrificing Specificity. RT² SYBR Green qPCR Master Mixes and Competitor I Master Mixes were used in qPCR assays to detect the human MMP13 and MMP15 mRNA in reference RNA. RT² SYBR Green qPCR Master mixes provide detection of the genes at an earlier threshold cycle value (C). The real-time dissociation curves and agarose gel electrophoresis characterization reveal the presence of a non-specific secondary product generated with the competitor's master mix which is not amplified by the RT² SYBR Green qPCR Master Mix.

RT² SYBR Green qPCR Master Mix from SABiosciences

- Proven in > 20,000 Genes & in > 600 Peer-Reviewed Publications

Compatible PCR Instruments	
Applied Biosystems (ABI):	5700, 7000, 7300, 7500, 7500 FAST, 7700, 7900HT, StepOnePlus (96- and 384-well blocks)
Bio-Rad:	CFX96, CFX384, iCycler, iQ5, MyiQ, Chromo4, Opticon, Opticon 2
Stratagene:	Mx3000P, Mx3005P, Mx4000
Roche:	LightCycler 480 (96- and 384-well blocks)
Eppendorf:	Mastercycler ep realplex
TaKaRa:	TP800
Fluidigm:	BioMark

RT ² qPCR Master Mixes		
Product	Size	Catalog #
RT ² SYBR Green w/ ROX qPCR Master Mix <i>For All ABI and Stratagene Instruments</i>	2 Arrays	PA-012
	12 Arrays	PA-012-12
	24 Arrays	PA-012-24
	4 Arrays*	PA-012-8
	25 mL	PA-112
RT ² SYBR Green w/ Fluorescein qPCR Master Mix <i>For Bio-Rad Instruments</i>	2 Arrays	PA-011
	12 Arrays	PA-011-12
	24 Arrays	PA-011-24
	4 Arrays*	PA-011-8
	25 mL	PA-111
RT ² SYBR Green qPCR Master Mix <i>Without Reference Dye</i>	2 Arrays	PA-010
	12 Arrays	PA-010-12
	24 Arrays	PA-010-24
	4 Arrays*	PA-010-8
	25 mL	PA-110

* Plate format is 384-well.

FREE Web-Based PCR ARRAY Data Analysis Software

This integrated web-based software package for the PCR Array System automatically performs all $\Delta\Delta C_t$ based fold-change calculations from your uploaded raw threshold cycle data. Simply providing the array's catalog number annotates the results to the correct gene list. The web portal delivers results not only in a tabular format but also in scatter, volcano, cluster-gram, and multi-group plots. Perform any pair-wise comparison between groups of experimental replicates by defining your own fold-change and statistical significance thresholds, or compare all of the groups side-by-side. The web portal also helps you correctly interpret the genomic DNA, reverse transcription efficiency, and positive PCR control well data. Make your pathway-focused gene expression analysis quick and painless with the PCR Array System and the PCR Array Data Analysis Suite.

- **Simple:** Just upload your data and define your parameters*
- **Convenient:** No downloading or installation required
- **Publication-Ready Output:** Export all results as FREE EXCEL files or PNG image files

* EXCEL-based data analysis templates are available from our website.

INSTRUCTIONS

- 1 Upload your data in a simple EXCEL file format.
- 2 Define your housekeeping genes and experimental groups.
- 3 Choose an automatically generated data analysis result

Take a test run with pre-loaded sample data set today:

<http://www.SABiosciences.com/pcrarraydataanalysis.php>

OR

Join our next live webinar entitled: "PCR Array Data Analysis Tutorial" at:

<http://www.SABiosciences.com/seminarlist.php>




Figure 1: The Scatter Plot Compares Gene Expression Levels Between Two Experimental Conditions. The graph plots the \log_{10} of normalized gene expression levels in a control condition (x-axis) versus an experimental condition (y-axis). Symbols outside the gray area indicate fold-differences larger than a threshold that you can define. The red symbols in the upper-left corner readily identify up-regulated genes, and the green symbols in the lower right corner readily identify down-regulated genes.




Figure 2: The Volcano Plot Indicates the Statistical Significance of Gene Expression Changes. The x-axis plots the \log_{10} of the fold-differences, while the y-axis plots their p-values based on student's t-test of your replicate raw C_t data. The red and green symbols outside the gray area conveniently have the same meaning as the Scatter Plot. Symbols in the Volcano Plot above the blue line readily identify fold-differences at least as statistically significant as a threshold that you can define.

Service Core for PCR ARRAY Gene Expression Analysis

SABiosciences provides comprehensive PCR Array Gene Expression Analysis Services for all of your real-time PCR-based needs. SABiosciences' Service Team has years of experience in gene expression profiling with quantitative SYBR[®] Green real-time PCR, as well as RNA isolation of the quantity and quality required by real-time PCR experiments.

Why Use Our Gene Expression Analysis Services?

- No access to a real-time PCR instrument
- Lack the time or manpower to perform the microarray or PCR experiments needed to complete your project
- Prefer to run an optimized pilot project before bringing the technology in-house for a large scale project

Let our in-house experts serve as your "external core facility" to save you time, money, and effort in performing your expression analysis projects. The services are flexible and can be tailored to exactly meet your specific needs. Simply submit experimental samples and receive gene expression profiling results in a matter of days.

Features of the PCR ARRAY Service

- Free consultation on your experimental design to help you take advantage of our expertise
- Confirmation of RNA sample concentration and quality
- Unbiased reverse transcription of RNA into cDNA template, including elimination of genomic DNA contamination
- Detailed and flexible data analysis according to your needs
- Quick results
- Competitive and affordable prices

RNA Isolation Services

- **Standard:** Preparation from Cells, Tissue, and Blood
- **Specialized:** Preparation from FFPE, LCM, FNAB and Other Samples





References Listed by Pathway-Focused PCR Arrays

ANGIOGENESIS

Airoldi I, Di Carlo E, Cocco C, Taverniti G, D'Antuono T, Ognio E, Watanabe M, Ribatti D, Pistoia V. Endogenous IL-12 triggers an antiangiogenic program in melanoma cells. *Proc Natl Acad Sci U S A*. 2007 Mar 6;104(10):3996-4001.

Basile DP, Fredrich K, Chelladurai B, Leonard EC, Parrish AR. Renal ischemia reperfusion inhibits VEGF expression and induces ADAMTS-1, a novel VEGF inhibitor. *Am J Physiol Renal Physiol*. 2008 Apr;294(4):F928-36.

Rangel-Moreno J, Hartson L, Navarro C, Gaxiola M, Selman M, Randall TD. Inducible bronchus-associated lymphoid tissue (iBALT) in patients with pulmonary complications of rheumatoid arthritis. *J Clin Invest*. 2006 Dec;116(12):3183-94.

van Beem RT, Noort WA, Voermans C, Kleijer M, ten Brinke A, van Ham SM, van der Schoot CE, Zwaginga JJ. The presence of activated CD4(+) T cells is essential for the formation of colony-forming unit-endothelial cells by CD14(+) cells. *J Immunol*. 2008 Apr 1;180(7):5141-8.

Brant KA, Fabisiak JP. Nickel and the microbial toxin, MALP-2, stimulate proangiogenic mediators from human lung fibroblasts via a HIF-1alpha and COX-2-mediated pathway. *Toxicol Sci*. 2009 Jan;107(1):227-37.

Liu Z, Kobayashi K, van Dinther M, van Heiningen SH, Valdimarsdottir G, van Laar T, Scharpfenecker M, Lowik CW, Goumans MJ, Ten Dijke P, Pardali E. VEGF and inhibitors of TGFbeta type-I receptor kinase synergistically promote J Cell Sci. 2009 Sep 15;122(Pt 18):3294-302.

McElroy MK, Kaushal S, Tran Cao HS, Moossa AR, Talamini MA, Hoffman RM, Bouvet M. Upregulation of thrombospondin-1 and angiogenesis in an aggressive human Mol Cancer Ther. 2009 Jul;8(7):1779-86.

APOPTOSIS

Dickson EW, Hogrefe CP, Ludwig PS, Ackermann LW, Stoll LL, Denning GM. Exercise Enhances Myocardial Ischemic Tolerance via an Opioid Receptor-Dependent Mechanism. *Am J Physiol Heart Circ Physiol*. 2007 Oct 19.

Lai YM, Mohammed KA, Nasreen N, Baumuratov A, Bellew BF, Antony VB. Induction of cell cycle arrest and apoptosis by BCG infection in cultured human bronchial airway epithelial cells. *Am J Physiol Lung Cell Mol Physiol*. 2007 Aug;293(2):L393-401.

Raeder S, Utheim TP, Utheim OA, Nicolaisen B, Roald B, Cai Y, Haug K, Kvalheim A, Messelt EB, Drolsum L, Reed JC, Lyberg T. Effects of organ culture and Optisol-GS storage on structural integrity, phenotypes, and apoptosis in cultured corneal epithelium. *Invest Ophthalmol Vis Sci*. 2007 Dec;48(12):5484-93.

Zhan Q, Signoretti S, Whitaker-Menezes D, Friedman TM, Korngold R, Murphy GF. Cytokeratin15-positive basal epithelial cells targeted in graft-versus-host disease express a constitutive antiapoptotic phenotype. *J Invest Dermatol*. 2007 Jan;127(1):106-15.

Mayoral R, Mollá B, Flores JM, Boscá L, Casado M, Martín Sanz P. Constitutive expression of cyclooxygenase 2 transgene in hepatocytes protects against liver injury. *Biochem J*. 2008 Aug 1.

McEuen JG, Beck SG, Bale TL. Failure to mount adaptive responses to stress results in dysregulation and cell death in the midbrain raphe. *J Neurosci*. 2008 Aug 13;28(33):8169-77.

Bose RN, Maurmann L, Mishur RJ, Yasui L, Gupta S, Grayburn WS, Hofstetter H, Milton T. Non-DNA-binding platinum anticancer agents: Cytotoxic activities of platinum-phosphato complexes towards human ovarian cancer cells. *Proc Natl Acad Sci U S A*. 2008 Nov 25;105(47):18314-9.

Uetani T, Nakayama H, Okayama H, Okura T, Higaki J, Inoue H, Higashiyama S. Insufficiency of proHB-EGF shedding enhances hypoxic cell death in H9c2 cardiomyoblasts via the activation of Caspase-3 and c-JUN N-terminal kinase. *J Biol Chem*. 2009 Feb 4.

Bok K, Prikhodko VG, Green KY, Sosnovtsev SV. Apoptosis in Murine Norovirus Infected RAW264.7 Cells is Associated with Survivin Downregulation. *J Virol*. 2009 Feb 11.

Mhyre AJ, Marcondes AM, Spaulding EY, Deeg HJ. Stroma-dependent apoptosis in clonal hematopoietic precursors correlates with expression of PYCARD. *Blood*. 2009 Jan 15;113(3):649-58.

Pru JK, Kaneko-Tarui T, Jurisicova A, Kashiwagi A, Selesniemi K, Tilly JL. Induction of proapoptotic gene expression and recruitment of p53 herald ovarian Reprod Sci. 2009 Apr;16(4):347-56.

ATHEROSCLEROSIS

Manning-Tobin JJ, Moore KJ, Seimon TA, Bell SA, Sharuk M, Alvarez-Leite JI, de Winther MP, Tabas I, Freeman MW. Loss of SR-A and CD36 Activity Reduces Atherosclerotic Lesion Complexity Without Abrogating Foam Cell Formation in Hyperlipidemic Mice. *Arterioscler Thromb Vasc Biol*. 2008 Oct 23.

BREAST CANCER & ESTROGEN RECEPTOR SIGNALING

Cooper C, Guo J, Yan Y, Chooniedass-Kothari S, Hube F, Hamedani MK, Murphy LC, Myal Y, Leygue E. Increasing the relative expression of endogenous non-coding Steroid Receptor RNA. *Nucleic Acids Res*. 2009 Jul;37(13):4518-31.

Adams BD, Cowee DM, White BA. The role of miR-206 in the epidermal growth factor (EGF) induced repression of Mol Endocrinol. 2009 Aug;23(8):1215-30.

APPLICATION: CANCER BIOLOGY




Figure: PIK3CA mRNA Expression in Multiple Lung Cancer Cell Lines. PIK3CA mRNA expression was compared among cell lines having different features, such as PIK3CA alterations or mutations of other genes involved in the EGFR signaling pathway. PIK3CA mRNA expression levels were expressed relative to the mean levels in six HBEC cell lines. PIK3CA mRNA expression in PIK3CA gain or EGFR mutant cell lines were significantly increased compared with that of wild-type cell lines. However, PIK3CA mutant lines do not express increased mRNA levels. Horizontal bars indicate mean values. The Kruskal-Wallis test with Dunn's multiple comparison test was used to determine significance.

Yamamoto, H., et al. *Cancer Research* 2008; 68: 6913-6921.

CANCER PATHWAYFINDER™

Gridley DS, Slater JM, Luo-Owen X, Rizvi A, Chapes SK, Stodieck LS, Ferguson VL, Pecaut MJ. Spaceflight effects on T lymphocyte distribution, function and gene expression. *J Appl Physiol*. 2009 Jan;106(1):194-202.

Das KK, Bajpai M, Kong Y, Liu J, Geng X, Das KM. Mesalamine suppresses the expression of TC22, a novel tropomyosin isoform *Mol Pharmacol*. 2009 Jul;76(1):183-91.

CELL CYCLE

Chen N, Ye XC, Chu K, Navone NM, Sage EH, Yu-Lee LY, Logothetis CJ, Lin SH. A secreted isoform of ErbB3 promotes osteonectin expression in bone and enhances the invasiveness of prostate cancer cells. *Cancer Res*. 2007 Jul 15;67(14):6544-8.

Chen S, Sims GP, Chen XX, Gu YY, Chen S, Lipsky PE. Modulatory effects of 1,25-dihydroxyvitamin D3 on human B cell differentiation. *J Immunol*. 2007 Aug 1;179(3):1634-47.

Lu SY, Sontag DP, Detillieux KA, Cattini PA. FGF-16 is released from neonatal cardiac myocytes and alters growth-related signaling: a possible role in postnatal development. *Am J Physiol Cell Physiol*. 2008 May;294(5):C1242-9.

Higgins S, Wong SH, Richner M, Rowe CL, Newgreen DF, Werther GA, Russo VC. Fibroblast growth factor 2 reactivates G1 checkpoint in SK-N-MC cells via *Endocrinology*. 2009 Sep;150(9):4044-55.

CHEMOKINES & RECEPTORS

Klopp AH, Spaeth EL, Dembinski JL, Woodward WA, Munshi A, Meyn RE, Cox JD, Andreoff M, Marini FC. Tumor irradiation increases the recruitment of circulating mesenchymal stem cells into the tumor microenvironment. *Cancer Res*. 2007 Dec 15;67(24):11687-95.

Miselis NR, Wu ZJ, Van Rooijen N, Kane AB. Targeting tumor-associated macrophages in an orthotopic murine model of diffuse malignant mesothelioma. *Mol Cancer Ther*. 2008 Apr;7(4):788-99.

Sinha S, Subramanian S, Proctor TM, Kaler LJ, Grafe M, Dahan R, Huan J, Vandenbark AA, Burrows GG, Offner H. A promising therapeutic approach for multiple sclerosis: recombinant T-cell receptor ligands modulate experimental autoimmune encephalomyelitis by reducing interleukin-17 production and inhibiting migration of encephalitogenic cells into the CNS *J Neurosci*. 2007 Nov 14;27(46):12531-9.

Sze SK, de Kleijn DP, Lai RC, Khia Way Tan E, Zhao H, Yeo KS, Low TY, Lian Q, Lee CN, Mitchell W, El Oakley RM, Lim SK. Elucidating the secretion proteome of human embryonic stem cell-derived mesenchymal stem cells. *Mol Cell Proteomics*. 2007 Oct;6(10):1680-9.

Lee PY, Kumagai Y, Li Y, Takeuchi O, Yoshida H, Weinstein J, Kellner ES, Nacionales D, Barker T, Kelly-Scumpia K, van Rooijen N, Kumar H, Kawai T, Satoh M, Akira S, Reeves WH. TLR7-dependent and FcγR-independent production of type I interferon in experimental mouse lupus. *J Exp Med*. 2008 Dec 22;205(13):2995-3006.

Lim J, Derrick SC, Kolibak K, Yang AL, Porcelli S, Jacobs WR, Morris SL. Early pulmonary cytokine and chemokine responses in mice immunized with three different vaccines against *Mycobacterium tuberculosis* determined by PCR array. *Clin Vaccine Immunol*. 2009 Jan;16(1):122-6.

Fukami N, Ramachandran S, Saini D, Walter M, Chapman W, Patterson GA, Mohanakumar T. Antibodies to MHC class I induce autoimmunity: role in the pathogenesis of chronic rejection. *J Immunol*. 2009 Jan 1;182(1):309-18.

Cheung KP, Yang E, Goldrath AW. Memory-like CD8+ T cells generated during homeostatic proliferation defer to *J Immunol*. 2009 Sep 1;183(5):3364-72.

Sundararaj KP, Samuvel DJ, Li Y, Sanders JJ, Lopes-Virella MF, Huang Y. Interleukin-6 released from fibroblasts is essential for up-regulation of matrix *J Biol Chem*. 2009 May 15;284(20):13714-24.

COMMON CYTOKINES

Gude N, Muraski J, Rubio M, Kajstura J, Schaefer E, Anversa P, Sussman MA. Akt promotes increased cardiomyocyte cycling and expansion of the cardiac progenitor cell population. *Circ Res*. 2006 Aug 18;99(4):381-8.

Keskin DB, Allan DS, Rybalov B, Andzelm MM, Stern JN, Kopcow HD, Koopman LA, Strominger JL. TGFβ promotes conversion of CD16+ peripheral blood NK cells into CD16- NK cells with similarities to decidual NK cells. *Proc Natl Acad Sci U S A*. 2007 Feb 27;104(9):3378-83.

Silver RF, Walrath J, Lee H, Jacobson BA, Horton H, Bowman MR, Nocka K, Sypek JP. Human Alveolar Macrophage Gene Responses to *Mycobacterium tuberculosis* Strains H37Ra and H37Rv. *Am J Respir Cell Mol Biol*. 2008 Sep 11.

Campeau PM, Rafei M, Boivin MN, Sun Y, Grabowski GA, Galipeau J. Characterization of Gaucher disease bone marrow mesenchymal stromal cells reveals *Blood*. 2009 Oct 8;114(15):3181-90.

Myskiw C, Arsenio J, van Bruggen R, Deschambault Y, Cao J. Vaccinia virus E3 suppresses expression of diverse cytokines through inhibition *J Virol*. 2009 Jul;83(13):6757-68.

CUSTOM PCR ARRAYS

Bowers RR, Kim JW, Otto TC, Lane MD. Stable stem cell commitment to the adipocyte lineage by inhibition of DNA methylation: Role of the BMP-4 gene. *Proc Natl Acad Sci U S A*. 2006 Aug 17.

Torcia MG, Santarlasci V, Cosmi L, Clemente A, Maggi L, Mangano VD, Verra F, Bancone G, Nebie I, Sirima BS, Liotta F, Frosali F, Angeli R, Severini C, Sannella AR, Bonini P, Lucibello M, Maggi E, Garaci E, Coluzzi M, Cozzolino F, Annunziato F, Romagnani S, Modiano D. Functional deficit of T regulatory cells in Fulani, an ethnic group with low susceptibility to *Plasmodium falciparum* malaria. *Proc Natl Acad Sci U S A*. 2008 Jan 15;105(2):646-51.

Shivapurkar N, Stastny V, Okumura N, Girard L, Xie Y, Prinsen C, Thunnissen FB, Wistuba II, Czerniak B, Frenkel E, Roth JA, Liloglou T, Xinarianos G, Field JK, Minna JD, Gazdar AF. Cytochrome b5, the newest member of the globin family, functions as a tumor suppressor gene. *Cancer Res*. 2008 Sep 15;68(18):7448-56.

Bousquet L, Pruvost A, Guyot AC, Farinotti R, Mabondzo A. Combination of Tenofovir and Emtricitabine plus Efavirenz: In Vitro Modulation of ABC transporter and Intracellular Drug Accumulation. *Antimicrob Agents Chemother*. 2008 Dec 15.

Huang J, Chen K, Huang J, Gong W, Dunlop NM, Howard OM, Bian X, Gao Y, Wang JM. Regulation of the leucocyte chemoattractant receptor FPR in glioblastoma cells by cell differentiation. *Carcinogenesis*. 2009 Feb;30(2):348-55.

Chahrour M, Jung SY, Shaw C, Zhou X, Wong ST, Qin J, Zoghbi HY. MeCP2, a Key Contributor to Neurological Disease, Activates and Represses Transcription Science. 2008 May 30;320(5880):1224-9.

Yamamoto H, Shigematsu H, Nomura M, Lockwood WW, Sato M, Okumura N, Soh J, Suzuki M, Wistuba II, Fong KM, Lee H, Toyooka S, Date H, Lam WL, Minna JD, Gazdar AF. PIK3CA mutations and copy number gains in human lung cancers. *Cancer Res*. 2008 Sep 1;68(17):6913-21.

Birrell MA, De Alba J, Catley MC, Hardaker E, Wong S, Collins M, Clarke DL, Farrow SN, Willson TM, Collins JL, Belvisi MG. Liver X receptor agonists increase airway reactivity in a model of asthma via increasing airway smooth muscle growth. *J Immunol*. 2008 Sep 15;181(6):4265-71.

Peduto L, Dulauroy S, Lochner M, Spath GF, Morales MA, Cumano A, Eberl G. Inflammation recapitulates the ontogeny of lymphoid stromal cells. *J Immunol*. 2009 May 1;182(9):5789-99.

Ianzini F, Kosmacek EA, Nelson ES, Napoli E, Erenpreisa J, Kalejs M, Mackey MA. Activation of meiosis-specific genes is associated with depolyploidization of Cancer Res. 2009 Mar 15;69(6):2296-304.

Schreiber TH, Deyev VV, Rosenblatt JD, Podack ER. Tumor-induced suppression of CTL expansion and subjugation by gp96-Ig Cancer Res. 2009 Mar 1;69(5):2026-33.

DiMeo TA, Anderson K, Phadke P, Fan C, Perou CM, Naber S, Kuperwasser C. A novel lung metastasis signature links Wnt signaling with cancer cell Cancer Res. 2009 Jul 1;69(13):5364-73.

DNA DAMAGE

Hoffman AE, Zheng T, Ba Y, Zhu Y. The circadian gene NPAS2, a putative tumor suppressor, is involved in DNA damage response. *Mol Cancer Res*. 2008 Sep;6(9):1461-8.

Chesnokova V, Wong C, Zonis S, Gruszka A, Wawrowsky K, Ren SG, Benschlomo A, Yu R. Diminished pancreatic {beta} cell mass in securin-null mice is caused by {beta} cell apoptosis and senescence. *Endocrinology*. 2009 Feb 12.

Santisteban M, Reiman JM, Asiedu MK, Behrens MD, Nassar A, Kalli KR, Haluska P, Ingle JN, Hartmann LC, Manjili MH, Radisky DC, Ferrone S, Knutson KL. Immune-induced epithelial to mesenchymal transition in vivo generates breast Cancer Res. 2009 Apr 1;69(7):2887-95.

DRUG METABOLISM: PHASE I ENZYMES

Ning B, Dial S, Sun Y, Wang J, Yang J, Guo L Systematic and simultaneous gene profiling of 84 drug-metabolizing genes in primary human hepatocytes *J Biomol Screen*. 2008 Mar;13(3):194-201.

EXTRACELLULAR MATRIX & ADHESION MOLECULES

Collins M, Hurford R, Chao Y, Patrick C, Langford TD. Variability in Virus and in Host During HIV-1 Infection. *Horn of Africa Journal of AIDS*. 2006 Volume III, Issue 1:38-46.

Gelain F, Bottai D, Vescovi A, Zhang S. Designer self-assembling Peptide nanofiber scaffolds for adult mouse neural stem cell 3-dimensional cultures. *PLoS ONE*. 2006 Dec 27;1:e119.

Johansen LD, Naumanen T, Knudsen A, Westerlund N, Gromova I, Juntilla M, Nielsen C, Bottzauw T, Tolkovsky A, Westermarck J, Coffey ET, Jaattela M, Kallunki T. IKAP localizes to membrane ruffles with filamin A and regulates actin cytoskeleton organization and cell migration. *J Cell Sci*. 2008 Mar 15;121(Pt 6):854-64.

Nacu N, Luzina IG, Highsmith K, Lockatell V, Pochetuhon K, Cooper ZA, Gillmeister MP, Todd NW, Atamas SP Macrophages produce TGF-beta-induced (beta-ig-h3) following ingestion of apoptotic cells and regulate MMP14 levels and collagen turnover in fibroblasts *J Immunol*. 2008 Apr 1;180(7):5036-44.

Oumouna-Benachour K, Hans CP, Suzuki Y, Naura A, Datta R, Belmadani S, Fallon K, Woods C, Boulares AH. Poly(ADP-ribose) polymerase inhibition reduces atherosclerotic plaque size and promotes factors of plaque stability in apolipoprotein E-deficient mice: effects on macrophage recruitment, nuclear factor-kappaB nuclear translocation, and foam cell death. *Circulation*. 2007 May 8;115(18):2442-50.

Shinoda Y, Ogata N, Higashikawa A, Manabe I, Shindo T, Yamada T, Kugimiya F, Ikeda T, Kawamura N, Kawasaki Y, Tsushima K, Takeda N, Nagai R, Hoshi K, Nakamura K, Chung UI, Kawaguchi H. Kruppel-like factor 5 causes cartilage degradation through transactivation of matrix metalloproteinase 9. *J Biol Chem*. 2008 Jul 10.

Tapia A, Salamonsen LA, Manuelpillai U, Dimitriadis E. Leukemia inhibitory factor promotes human first trimester extravillous trophoblast adhesion to extracellular matrix and secretion of tissue inhibitor of metalloproteinases-1 and -2. *Hum Reprod*. 2008 Aug;23(8):1724-32.

Yen JH, Khayrullina T, Ganea D. PGE2-induced metalloproteinase-9 is essential for dendritic cell migration *Blood*. 2007 Oct 9.

Zauli G, Corallini F, Bossi F, Fischetti F, Durigutto P, Celeghini C, Tedesco F, Secchiero P. Osteoprotegerin increases leukocyte adhesion to endothelial cells both in vitro and in vivo. *Blood*. 2007 Jul 15;110(2):536-43.

Lucio-Eterovic AK, Piao Y, de Groot JF. Mediators of glioblastoma resistance and invasion during antivasculature endothelial *Clin Cancer Res*. 2009 Jul 15;15(14):4589-99.

Kuhn AR, Schlauch K, Lao R, Halayko AJ, Gerthoffer WT, Singer CA. MicroRNA Expression in Human Airway Smooth Muscle Cells: Role of miR-25 in Am *J Respir Cell Mol Biol*. 2009 Jun 18.

ENDOTHELIAL CELL BIOLOGY

Li Z, Suzuki Y, Huang M, Cao F, Xie X, Connolly AJ, Yang PC, Wu JC Comparison of Reporter Gene and Iron Particle Labeling for Tracking Fate of Human Embryonic Stem Cells and Differentiated Endothelial Cells in Living Subjects. *Stem Cells*. 2008 Apr;26(4):864-73.

Asada S, Takahashi T, Isodono K, Adachi A, Imoto H, Ogata T, Ueyama T, Matsubara H, Oh H. Downregulation of Dicer expression by serum withdrawal sensitizes human endothelial cells to apoptosis. *Am J Physiol Heart Circ Physiol*. 2008 Dec;295(6):H2512-21.

GROWTH FACTORS

Stirling DP, Liu S, Kubers P, Yong VW Depletion of Ly6G/Gr-1 leukocytes after spinal cord injury in mice alters wound healing and worsens neurological outcome. *J Neurosci*. 2009 Jan 21;29(3):753-64.

HYPOXIA SIGNALING

Wendler CC, Amatya S, McClaskey C, Ghatpande S, Fredholm BB, Rivkees SA. A1 adenosine receptors play an essential role in protecting the embryo against hypoxia. *Proc Natl Acad Sci U S A*. 2007 Jun 5;104(23):9697-702.

Mueller BR, Bale TL. Sex-specific programming of offspring emotionality after stress early in pregnancy. *J Neurosci*. 2008 Sep 3;28(36):9055-65.

APPLICATION: INFLAMMATORY CYTOKINES & RECEPTORS




Figure: Lineage-dependent Differences in Expression of Treg Functional Genes. Quantitative PCR analysis of selected genes related to Treg function in CD4⁺CD25⁺ cells from Mossi, Fulani, and European donors infected with *Plasmodium falciparum* (malaria). CD4⁺CD25⁺ cells were isolated from 12 Mossi (red bars) and 12 Fulani (blue bars) donors included in the study and 10 European donors (gray bars), and lysed to obtain total RNA. Equal amounts of RNA (50 ng) from each donor were reverse-transcribed and amplified in duplicate in RT² Custom PCR Arrays to simultaneously examine the mRNA levels of nine selected genes related to Treg activity using PPIA, GAPDH, and ACTB as housekeeping genes (HKG). Data was normalized to the mean values of the HKG, and a relative amount of RNA was calculated using the 2^{-ΔΔC} method.

Torcia MG, Santarasci V, et al. Proc Natl Acad Sci U S A. 2008 Jan 15;105:646-651.

Aubin E, Lemieux R, Bazin R. Absence of cytokine modulation following therapeutic infusion of intravenous immunoglobulin or anti-red blood cell antibodies in a mouse model of immune thrombocytopenic purpura. *Br J Haematol*. 2007 Mar;136(6):837-43.

Lakshmanan U, Porter AG. Caspase-4 interacts with TNF receptor-associated factor 6 and mediates lipopolysaccharide-induced NF-kappaB-dependent production of IL-8 and CC chemokine ligand 4 (macrophage-inflammatory protein-1). *J Immunol.* 2007 Dec 15;179(12):8480-90.

Lee JK, McCoy MK, Harms AS, Ruhn KA, Gold SJ, Tansey MG. Regulator of G-protein signaling 10 promotes dopaminergic neuron survival via regulation of the microglial inflammatory response. *J Neurosci.* 2008 Aug 20;28(34):8517-28.

Chittur S, Parr B, Marcovici G. Inhibition of Inflammatory Gene Expression in Keratinocytes Using a Composition Evid Based Complement Alternat Med. 2009 Aug 19.

Kawaguchi A, Orba Y, Kimura T, Iha H, Ogata M, Tsuji T, Aina A, Sata T, Okamoto T, Hall WW, Sawa H, Hasegawa H. Inhibition of the SDF-1alpha-CXCR4 axis by the CXCR4 antagonist AMD3100. *Blood.* 2009 Oct 1;114(14):2961-8.

Kanda T, Steele R, Ray R, Ray RB. Inhibition of intrahepatic gamma interferon production by hepatitis C virus. *J Virol.* 2009 Sep;83(17):8463-9.

Huston JM, Rosas-Ballina M, Xue X, Dowling O, Ochani K, Ochani M, Yeboah MM, Chatterjee PK, Tracey KJ, Metz CN. Cholinergic neural signals to the spleen down-regulate leukocyte trafficking via J Immunol. 2009 Jul 1;183(1):552-9.

MacNamara KC, Racine R, Chatterjee M, Borjesson D, Winslow GM. Diminished hematopoietic activity associated with alterations in innate and Infect Immun. 2009 Sep;77(9):4061-9.

Montgomery CP, Daum RS. Transcription of inflammatory genes in the lung after infection with Infect Immun. 2009 May;77(5):2159-67.

Choi K, Ahn YH, Gibbons DL, Tran HT, Creighton CJ, Girard L, Minna JD, Qin FX, Kurie JM. Distinct biological roles for the notch ligands Jagged-1 and Jagged-2. *J Biol Chem.* 2009 Jun 26;284(26):17766-74.

INFLAMMATORY RESPONSE & AUTOIMMUNITY

Rafei M, Campeau PM, Wu JH, Birman E, Forner K, Boivin MN, Galipeau J. Selective inhibition of CCR2 expressing lymphomyeloid cells in experimental J Immunol. 2009 Mar 1;182(5):2620-7.

INSULIN SIGNALING

Narvaez CJ, Matthews D, Broun E, Chan M, Welsh J. Lean phenotype and resistance to diet-induced obesity in vitamin D receptor knockout mice correlates with induction of uncoupling protein-1 in white adipose tissue. *Endocrinology.* 2009 Feb;150(2):651-61.

INTERFERON α , β RESPONSE

Tilton JC, Manion MM, Luskin MR, Johnson AJ, Patamawenu AA, Hallahan CW, Cogliano-Shutta NA, Mican JM, Davey RT Jr, Kottlilil S, Lifson JD, Metcalf JA, Lempicki RA, Connors M. Human immunodeficiency virus viremia induces plasmacytoid dendritic cell activation in vivo and diminished alpha interferon production in vitro. *J Virol.* 2008 Apr;82(8):3997-4006.

Li Q, Xu B, Michie SA, Rubins KH, Schreiber RD, McDevitt HO. Interferon-alpha initiates type 1 diabetes in nonobese diabetic mice. *Proc Natl Acad Sci U S A.* 2008 Aug 26;105(34):12439-44.

Christensen JE, Simonsen S, Fenger C, Sørensen MR, Moos T, Christensen JP, Finsen B, Thomsen AR. Fulminant lymphocytic choriomeningitis virus-induced inflammation of the CNS involves a cytokine-chemokine-cytokine-chemokine cascade. *J Immunol.* 2009 Jan 15;182(2):1079-87.

INTERFERON (IFN) & RECEPTOR

Yoshimi R, Chang TH, Wang H, Atsumi T, Morse HC 3rd, Ozato K. Gene disruption study reveals a nonredundant role for TRIM21/Ro52 in J Immunol. 2009 Jun 15;182(12):7527-38.

JAK / STAT SIGNALING

Albanesi C, Fairchild HR, Madonna S, Scarponi C, De Pit'0, Leung DY, Howell MD. IL-4 and IL-13 negatively regulate TNF-alpha- and IFN-gamma-induced beta-defensin expression through STAT-6, suppressor of cytokine signaling (SOCS)-1, and SOCS-3. *J Immunol.* 2007 Jul 15;179(2):984-92.

Duechting A, Tschöpe C, Kaiser H, Lamkemeyer T, Tanaka N, Aberle S, Lang F, Torresi J, Kandolf R, Bock CT. Human Parvovirus B19 NS1 Protein Modulates Inflammatory Signaling by Activation of STAT3/PIAS3 in Human Endothelial Cells. *J Virol.* 2008 Jun 11.

Noronha AM, Liang Y, Hetzel JT, Hasturk H, Kantarci A, Stucchi A, Zhang Y, Nikolajczyk BS, Farraye FA, Ganley-Leal LM. Hyperactivated B cells in human inflammatory bowel disease. *J Leukoc Biol.* 2009 Oct;86(4):1007-16.

LIPOPROTEIN SIGNALING & CHOLESTEROL

Zhou C, King N, Chen KY, Breslow JL. Activation of pregnane X receptor induces hypercholesterolemia in wild-type and J Lipid Res. 2009 May 12.

MAP KINASE SIGNALING

Ruetgger A, Schueler S, Mollenhauer JA, Wiederanders B. Cathepsins B, K and L are regulated by a defined collagen type II peptide via activation of classical protein kinase C and p38 MAP kinase in articular chondrocytes. *J Biol Chem.* 2007 Nov 8.

Singh S, Yin X, Pisano MM, Greene RM. Molecular profiles of mitogen activated protein kinase signaling pathways in orofacial development. *Birth Defects Res A Clin Mol Teratol.* 2007 Jan;79(1):35-44.

Makarenkova HP, Hoffman MP, Beenken A, Eliseenkova AV, Meech R, Tsau C, Patel VN, Lang RA, Mohammadi M. Differential interactions of FGFs with heparan sulfate control gradient formation. *Sci Signal.* 2009 Sep 15;2(88):ra55.

APPLICATION: NEUROSCIENCE




Figure: Gene Expression Changes in Hypothalamus of MECP2 Mouse Models. Validation of expression changes for 66 genes by qPCR analysis. Gene expression levels from microarray analyses were validated in four MECP2-Tg males and four MECP2-null males. Data is plotted as relative up-regulation (red) or down-regulation (blue) over WT ($P < 0.05$, t test). Each row represents a single gene, and each column represents data for four samples from each genotype.

Chahrour, M., et al. *Science* 2008; 320: 1224-1229.

NF κ B SIGNALING

Han YM, Romero R, Kim JS, Tarca AL, Kim SK, Draghici S, Kusanovic JP, Gotsch F, Mittal P, Hassan S, Kim CJ. Region-Specific Gene Expression Profiling: Novel Evidence for Biological Heterogeneity of the Human Amnion. *Biol Reprod.* 2008 Aug 6.

Singh K, Sinha S, Malonia SK, Bist P, Tergaonkar V, Chattopadhyay S. Tumor suppressor SMAR1 represses IkappaBalpha expression and inhibits p65 transactivation through matrix attachment regions. *J Biol Chem.* 2009 Jan 9;284(2):1267-78.

Farhana L, Dawson MI, Xu L, Fontana JA. SHP and Sin3A expression are essential for adamantyl-substituted retinoid-related Mol Cancer Ther. 2009 Jun;8(6):1625-35.

OXIDATIVE STRESS & ANTIOXIDANT DEFENSE

Chua PJ, Yip GW, Bay BH. Cell cycle arrest induced by hydrogen peroxide is associated with modulation of Exp Biol Med (Maywood). 2009 Sep;234(9):1086-94.

Xu Y, Saegusa C, Schehr A, Grant S, Whitsett JA, Ikegami M. C/EBP{alpha} is required for pulmonary cytoprotection during hyperoxia. Am J Physiol Lung Cell Mol Physiol. 2009 Aug;297(2):L286-98.

Baqai FP, Gridley DS, Slater JM, Luo-Owen X, Stodieck LS, Ferguson V, Chapes SK, Pecaut MJ. Effects of spaceflight on innate immune function and antioxidant gene expression. J Appl Physiol. 2009 Jun;106(6):1935-42.

Chen A, Cuevas I, Kenny PA, Miyake H, Mace K, Ghajar C, Boudreau A, Bissell M, Boudreau N. Endothelial cell migration and vascular endothelial growth factor expression are Cancer Res. 2009 Aug 15;69(16):6721-9.

Polytarchou C, Pfau R, Hatzia Apostolou M, Tsihli PN The JmjC domain histone demethylase Ndy1 regulates redox homeostasis and protects cells from oxidative stress. Mol Cell Biol. 2008 Dec;28(24):7451-64.

OSTEOGENESIS

Gouttenoire J, Valcourt U, Bougault C, Aubert-Foucher E, Arnaud E, Giraud L, Mallein-Gerin F. Knockdown of the intraflagellar transport protein IFT46 stimulates selective gene expression in mouse chondrocytes and affects early development in zebrafish. J Biol Chem. 2007 Oct 19;282(42):30960-73.

Wescott DC, Pinkerton MN, Gaffey BJ, Beggs KT, Milne TJ, Meikle MC. Osteogenic gene expression by human periodontal ligament cells under cyclic tension. J Dent Res. 2007 Dec;86(12):1212-6.

Flowers S, Nagl NG Jr, Beck GR Jr, Moran E Antagonistic roles for BRM and BRG1 SWI/SNF complexes in differentiation. J Biol Chem. 2009 Jan 14.

Meikle MC Northcroft Memorial Lecture 2007. A century of progress: advances in orthodontics since the foundation of the British Society for the Study of Orthodontics. J Orthod. 2008 Sep;35(3):176-90.

p53 SIGNALING

Brockschmidt C, Hirner H, Huber N, Eismann T, Hillenbrand A, Giamas G, Radunsky B, Ammerpohl O, Bohm B, Henne-Bruns D, Kalthoff H, Leithauser F, Trauzold A, Knippschild U Anti-apoptotic and growth-stimulatory functions of CK1 delta and epsilon in ductal adenocarcinoma of the pancreas are inhibited by IC261 in vitro and in vivo. Gut. 2008 Jun;57(6):799-806.

Cutts BA, Sjogren AK, Andersson KM, Wahlstrom AM, Karlsson C, Swolin B, Bergo MO. NF1 deficiency co-operates with oncogenic K-RAS to induce acute myeloid leukemia Blood. 2009 Aug 26.

Barkic M, Crnomarkovic S, Grabusic K, Bogetic I, Panic L, Tamarut S, Cokaric M, Jeric I, Vidak S, Volarevic S. The p53 tumor suppressor causes congenital malformations in Rpl24-deficient mice Mol Cell Biol. 2009 May;29(10):2489-504.

SIGNAL TRANSDUCTION PATHWAYFINDER

Hager-Theodorides AL, Furmanski AL, Ross SE, Outram SV, Rowbotham NJ, Crompton T. The Gli3 transcription factor expressed in the thymus stroma controls thymocyte J Immunol. 2009 Sep 1;183(5):3023-32.

Bojanic DD, Tarr PT, Gale GD, Smith DJ, Bok D, Chen B, Nusinowitz S, Lovgren-Sandblom A, Bjorkhem I, Edwards PA. Differential expression and function of ABCG1 and ABCG4 during development and J Lipid Res. 2009 Jul 26.

Rowbotham NJ, Hager-Theodorides AL, Furmanski AL, Ross SE, Outram SV, Dessens JT, Crompton T. Sonic hedgehog negatively regulates pre-TCR-induced differentiation by a Blood. 2009 May 21;113(21):5144-56.

T CELL ANERGY & IMMUNE TOLERANCE

Gorczynski RM, Chen Z, He W, Khatri I, Sun Y, Yu K, Boudakov I. Expression of a CD200 transgene is necessary for induction but not maintenance of J Immunol. 2009 Aug 1;183(3):1560-8.

TGFβ BMP SIGNALING

Garamszegi N, Garamszegi SP, Shehadeh LA, Scully SP. Extracellular matrix-induced gene expression in human breast cancer cells. Mol Cancer Res. 2009 Mar;7(3):319-29.

Th1, Th2, Th3

Jacob N, Yang H, Pricop L, Liu Y, Gao X, Zheng SG, Wang J, Gao HX, Putterman C, Koss MN, Stohl W, Jacob CO. Accelerated pathological and clinical nephritis in systemic lupus erythematosus-prone New Zealand Mixed 2328 mice doubly deficient in TNF receptor 1 and TNF receptor 2 via a Th17-associated pathway. J Immunol. 2009 Feb 15;182(4):2532-41.

Roman J, Rangasamy T, Guo J, Sugunan S, Meednu N, Packirisamy G, Shimoda L, Golding A, Semenza G, Georas SN. T Cell Activation Under Hypoxic Conditions Enhances Interferon-gamma Secretion. Am J Respir Cell Mol Biol. 2009 Apr 16.

Th17 FOR AUTOIMMUNITY & INFLAMMATION

Lochner M, Peduto L, Cherrier M, Sawa S, Langa F, Varona R, Riethmacher D, Si-Tahar M, Di Santo JP, Eberl G. In vivo equilibrium of proinflammatory IL-17+ and regulatory IL-10+ Foxp3+ RORgamma+ T cells. J Exp Med. 2008 Jun 9;205(6):1381-93.

TOLL-LIKE RECEPTOR SIGNALING

Brant KA, Fabisiak JP. Nickel Alterations of TLR2-Dependent Chemokine Profiles in Lung Fibroblasts are Mediated by COX-2. Am J Respir Cell Mol Biol. 2007 Dec 20.

Ornatowska M, Azim AC, Wang X, Christman JW, Xiao L, Joo M, Sadikot RT. Functional genomics of silencing TREM-1 on TLR4 signaling in macrophages. Am J Physiol Lung Cell Mol Physiol. 2007 Dec;293(6):L1377-84.

Tomchuck SL, Zvezdaryk KJ, Coffelt SB, Waterman RS, Danka ES, Scandurro AB Toll-Like Receptors on Human Mesenchymal Stem Cells Drive their Migration and Immunomodulating Responses. Stem Cells. 2007 Oct 4.

West J, Damania B Upregulation of the TLR3 pathway by Kaposi's sarcoma-associated herpesvirus during primary infection. J Virol. 2008 Jun;82(11):5440-9.

Prêle CM, Woodward EA, Bislely J, Keith-Magee A, Nicholson SE, Hart PH. SOCS1 regulates the IFN but not NFκB pathway in TLR-stimulated human monocytes and macrophages. J Immunol. 2008 Dec 1;181(11):8018-26.

Shin H, Zhang Y, Jagannathan M, Hasturk H, Kantarci A, Liu H, Van Dyke TE, Ganley-Leal LM, Nikolajczyk BS. B cells from periodontal disease patients express surface Toll-like receptor 4 J Leukoc Biol. 2008 Dec 31.

Samuel DJ, Sundararaj KP, Nareika A, Lopes-Virella MF, Huang Y. Lactate boosts TLR4 signaling and NF-κB pathway-mediated gene transcription in macrophages via monocarboxylate transporters and MD-2 up-regulation. J Immunol. 2009 Feb 15;182(4):2476-84.

Chen GY, Shiah HC, Su HJ, Chen CY, Chuang YJ, Lo WH, Huang JL, Chuang CK, Hwang SM, Hu YC. Baculovirus transduction of mesenchymal stem cells triggers the toll-like J Virol. 2009 Oct;83(20):10548-56.

Zhang X, Alnaeeli M, Singh B, Teng YT. Involvement of SOCS3 in regulation of CD11c+ dendritic cell-derived Infect Immun. 2009 May;77(5):2000-9.

TRANSCRIPTION FACTORS

Yuvaraj S, Griffin AC, Sundaram K, Kirkwood KL, Norris JS, Reddy SV. A novel function of CXCL13 to stimulate RANK ligand expression in oral squamous Mol Cancer Res. 2009 Aug;7(8):1399-407.

Middeldorp J, Kamphuis W, Sluijs JA, Achoui D, Leenaars CH, Feenstra MG, van Tijn P, Fischer DF, Berkens C, Ovaas H, Quinlan RA, Hol EM. Intermediate filament transcription in astrocytes is repressed by proteasome FASEB J. 2009 Aug;23(8):2710-26.

TOXICOLOGY

Joseph, P. et al. Heme-oxygenase 1 Gene Expression is a Marker for Hexavalent Chromium-Induced Stress and Toxicity in Human Dermal Fibroblasts. Toxicol. Sci., Jun 2008; 103: 325-334.

Marshall, N.B, et al. Functional Characterization and Gene Expression Analysis of CD4+CD25+ Regulatory T Cells Generated in Mice Treated with 2,3,7,8-Tetrachlorodibenzo-p-Dioxin. J. Immunol., Aug 2008; 181: 2382 - 2391.

TUMOR METASTASIS

Wang J, Shiozawa Y, Wang J, Wang Y, Jung Y, Pienta KJ, Mehra R, Lober R, Taichman RS. The role of CXCR7/RDC1 as a chemokine receptor for CXCL12/ SDF-1 in prostate cancer. J Biol Chem. 2007 Dec 5.

WNT SIGNALING

Aslanidi G, Kroutov V, Philipsberg G, Lamb K, Campbell-Thompson M, Walter GA, Kurenov S, Ignacio Aguirre J, Keller P, Hankenson K, Macdougald OA, Zolotukhin S. Ectopic expression of Wnt10b decreases adiposity and improves glucose homeostasis in obese rats. Am J Physiol Endocrinol Metab. 2007 Sep;293(3):E726-36.

Giuliani N, Morandi F, Tagliaferri S, Lazzaretti M, Donofrio G, Bonomini S, Sala R, Mangoni M, Rizzoli V. Production of wnt inhibitors by myeloma cells: potential effects on canonical wnt pathway in the bone microenvironment. Cancer Res. 2007 Aug 15;67(16):7665-74.

Guo W, Flanagan J, Jasuja R, Kirkland J, Jiang L, Bhasin S. The Effects of Myostatin on Adipogenic Differentiation of Human Bone Marrow-derived Mesenchymal Stem Cells Are Mediated through Cross-communication between Smad3 and Wnt/{beta}-Catenin Signaling Pathways. J Biol Chem. 2008 Apr 4;283(14):9136-45.

Hidestrand M, Richards-Malcolm S, Gurley CM, Nolen G, Grimes B, Waterstrat A, Zant GV, Peterson CA. Sca-1-expressing nonmyogenic cells contribute to fibrosis in aged skeletal muscle. J Gerontol A Biol Sci Med Sci. 2008 Jun;63(6):566-79.

Welters HJ, Oknianska A, Erdmann KS, Ryffel GU, Morgan NG. The protein tyrosine phosphatase-BL, modulates pancreatic beta-cell proliferation by interaction with the Wnt signaling pathway. J Endocrinol. 2008 Jun;197(3):543-52.

Malhotra S, Baba Y, Garrett KP, Staal FJ, Gerstein R, Kincade PW. Contrasting responses of lymphoid progenitors to canonical and noncanonical Wnt signals. J Immunol. 2008 Sep 15;181(6):3955-64.

Gelebart P, Anand M, Armanious H, Peters AC, Dien Bard J, Amin HM, Lai R. Constitutive activation of the Wnt canonical pathway in mantle cell lymphoma. Blood. 2008 Dec 15;112(13):5171-9.

Jonassen JA, San Agustin J, Follit JA, Pazour GJ. Deletion of IFT20 in the mouse kidney causes misorientation of the mitotic spindle and cystic kidney disease. J Cell Biol. 2008 Nov 3;183(3):377-84.

Wu QL, Zierold C, Ranheim EA. Dysregulation of frizzled 6 is a critical component of B cell leukemogenesis in a mouse model of chronic lymphocytic leukemia. Blood. 2009 Jan 28.

Ghosh MC, Collins GD, Vandanmagsar B, Patel K, Brill M, Carter A, Lustig A, Becker KG, Wood WW 3rd, Emeche CD, French AD, O'Connell MP, Xu M, Weeraratna AT, Taub DD. Activation of Wnt5A signaling is required for CXC chemokine ligand 12-mediated Blood. 2009 Aug 13;114(7):1366-73.

Hussain M, Rao M, Humphries AE, Hong JA, Liu F, Yang M, Caragacianu D, Schrumpp DS. Tobacco smoke induces polycomb-mediated repression of Dickkopf-1 in lung cancer Cancer Res. 2009 Apr 15;69(8):3570-8.

Data is used with permission from Science, PNAS, Cancer Research, & respective authors.

NOTICE TO PURCHASERS

This product is intended for research purposes only and is not intended for drug or diagnostic purposes or for human use. Use of kit components for reproduction of any primer pair mix, to modify kit components for resale or to use RT² Profiler PCR Arrays, RT² qPCR Primer Assays, RT² miRNA PCR Arrays & Assays, Methyl-Profiler PCR Array & qPCR Primer Assays, and ChampionChIP PCR Arrays & Real-Time PCR Primers to manufacture commercial products without written approval of SABiosciences Corporation is expressly prohibited.

PRODUCT WARRANTY

This warranty limits our liability to the replacement of this product in the event the product fails to perform due to any manufacturing defect. SABiosciences Corporation makes no other warranties of any kind, expressed or implied, including without limitation, warranties of merchantability or fitness for a particular purpose. SABiosciences Corporation shall not be liable for any direct, indirect, consequential or incidental damages arising out of the use, the results of use or the inability to use this product.

LIMITED LICENSE STATEMENTS

Use of this product is covered by one of more of the following US patents and corresponding patent claims outside the US: 5,079,352, 5,789,224, 5,618,711, 6,127,155, 5,677,152 (claims 1 to 23 only), 5,773,258 (claims 1 and 6 only), 5,407,800, 5,322,770, 5,310,652, 5,994,056, 6,171,785, and claims outside the US corresponding to US Patent No. 4,889,818. The purchase of this product includes a limited, non-transferable immunity from suit under the foregoing patent claims for using only this amount of product for the purchaser's own internal research. No right under any other patent claim (such as apparatus or system claims in US Patent No. 6,814,934) and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is conveyed expressly, by implication, or by estoppel. This product is for research use only. Diagnostic uses under Roche patents require a separate license from Roche. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

This product is provided under an agreement between Molecular Probes, Inc. and SABiosciences and the manufacture, use, sale, or import of this product is subject to one or more of U.S. Patent Nos. 5,436,134; 5,658,751 and corresponding international equivalents, owned by Invitrogen Corp. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer, where such research does not include testing, analysis or screening services for any third party in return for compensation on a per test basis. The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. For information on purchasing a license to this product for purposes other than research, contact Molecular Probes, Inc., Business Development, 29851 Willow Creek Road, Eugene, OR 97402 Tel: (541) 465-8300, Fax: (541) 335-0504.

The purchase of this product includes a limited, non-transferable license under specific claims of U.S. Patent Nos. 6,174,670, 6,569,627 and 5,871,908, owned by the University of Utah Research Foundation or Evotec Biosystems GmbH and licensed to Idaho Technology, Inc. and Roche Diagnostics GmbH, to use only the enclosed amount of product according to the specified protocols. No right is conveyed, expressly, by implication, or by estoppel, to use any instrument or system under any claim of U.S. Patent Nos. 6,174,670, 6,569,627 and 5,871,908, other than for the amount of product contained herein.

Patent Pending.

COMPLETE PCR ARRAY PRODUCT CATALOG

RT ² Profiler PCR ARRAYS*	Human	Mouse	Rat
Alzheimer's Disease	PAHS-057	PAMM-057	
Angiogenesis	PAHS-024	PAMM-024	PARN-024
Angiogenic Growth Factors and Inhibitors	PAHS-072	PAMM-072	PARN-072
Apoptosis	PAHS-012	PAMM-012	
Atherosclerosis	PAHS-038	PAMM-038	
Autophagy	PAHS-084	PAMM-084	
Breast Cancer and Estrogen Receptor Signaling	PAHS-005	PAMM-005	PARN-005
cAMP / Calcium Signaling PathwayFinder™	PAHS-066	PAMM-066	
Cancer Drug Resistance and Metabolism	PAHS-004	PAMM-004	
Cancer PathwayFinder™	PAHS-033	PAMM-033	PARN-033
Cell Cycle	PAHS-020	PAMM-020	PARN-020
Cell Surface Markers	PAHS-055		
Chemokines and Receptors	PAHS-022	PAMM-022	PARN-022
Common Cytokines	PAHS-021	PAMM-021	PARN-021
Dendritic and Antigen Presenting Cell	PAHS-406		
Diabetes	PAHS-023	PAMM-023	PARN-023
DNA Damage Signaling Pathway	PAHS-029	PAMM-029	PARN-029
Drug Metabolism	PAHS-002	PAMM-002	PARN-002
Drug Metabolism: Phase I Enzymes	PAHS-068		PARN-068
Drug Metabolism: Phase II Enzymes	PAHS-069		PARN-069
Drug Transporters	PAHS-070		PARN-070
EGF / PDGF Signaling Pathway	PAHS-040		PARN-040
Embryonic Stem Cells	PAHS-081	PAMM-081	
Endothelial Cell Biology	PAHS-015	PAMM-015	PARN-015
Extracellular Matrix and Adhesion Molecules	PAHS-013	PAMM-013	PARN-013
GPCR Signaling PathwayFinder™	PAHS-071		PARN-071
Growth Factors	PAHS-041	PAMM-041	PARN-041
Heat Shock Proteins	PAHS-076	PAMM-076	
Hedgehog Signaling Pathway	PAHS-078	PAMM-078	
Hematopoietic Stem Cells and Hematopoiesis	PAHS-054		
HIV Infection and Host Response	PAHS-051		
Homeobox (HOX) Genes	PAHS-083	PAMM-083	
Housekeeping Genes	PAHS-000	PAMM-000	PARN-000
Hypoxia Signaling Pathway	PAHS-032	PAMM-032	PARN-032
Inflammatory Cytokines and Receptors	PAHS-011	PAMM-011	PARN-011
Inflammatory Response and Autoimmunity	PAHS-077		
Innate and Adaptive Immune Responses	PAHS-052	PAMM-052	
Insulin Signaling Pathway	PAHS-030	PAMM-030	PARN-030
Interferon α,β Response	PAHS-016	PAMM-016	
Interferons and Receptors	PAHS-064	PAMM-064	
JAK / STAT Signaling Pathway	PAHS-039	PAMM-039	PARN-039
Lipoprotein Signaling and Cholesterol Metabolism	PAHS-080	PAMM-080	PARN-080
MAP Kinase Signaling Pathway	PAHS-061	PAMM-061	PARN-061
Mesenchymal Stem Cells	PAHS-082	PAMM-082	
Neurogenesis and Neural Stem Cell	PAHS-404	PAMM-404	
Neuroscience Ion Channels and Transporters	PAHS-036	PAMM-036	PARN-036
Neurotransmitter Receptors and Regulators	PAHS-060	PAMM-060	PARN-060
Neurotrophins and Receptors	PAHS-031	PAMM-031	PARN-031
NFκB Signaling Pathway	PAHS-025	PAMM-025	PARN-025
Nitric Oxide Signaling Pathway	PAHS-062	PAMM-062	PARN-062
Notch Signaling Pathway	PAHS-059	PAMM-059	
Nuclear Receptors and Coregulators	PAHS-056	PAMM-056	PARN-056
Obesity		PAMM-017	PARN-017
Osteogenesis	PAHS-026	PAMM-026	PARN-026
Oxidative Stress and Antioxidant Defense	PAHS-065	PAMM-065	PARN-065
p53 Signaling Pathway	PAHS-027	PAMM-027	PARN-027
PI3K-AKT Signaling Pathway	PAHS-058		
Signal Transduction PathwayFinder™	PAHS-014	PAMM-014	PARN-014
Stem Cell	PAHS-405	PAMM-405	PARN-405
Stress and Toxicity PathwayFinder™	PAHS-003	PAMM-003	PARN-003
Stress Response to Cellular Damage		PAMM-019	
T-Cell and B-Cell Activation	PAHS-053	PAMM-053	
T Cell Anergy and Immune Tolerance		PAMM-074	
TGFβ / BMP Signaling Pathway	PAHS-035	PAMM-035	PARN-035
Th1-Th2-Th3	PAHS-034	PAMM-034	PARN-034
Th17 for Autoimmunity and Inflammation	PAHS-073	PAMM-073	
TNF Ligands and Receptors	PAHS-063	PAMM-063	
Toll-Like Receptor Signaling Pathway	PAHS-018	PAMM-018	PARN-018
Transcription Factors	PAHS-075		
Tumor Metastasis	PAHS-028	PAMM-028	
Ubiquitination (Ubiquitylation)	PAHS-079	PAMM-079	
Wnt Signaling Pathway	PAHS-043	PAMM-043	PARN-043

PCR ARRAYS	Size
96-well	2 Arrays
96-well	12 Arrays
96-well	24 Arrays
384-well	4 Arrays
Add up to 4 Genes to any PCR Array	

Compatible Instrument	Plate
Applied Biosystems (ABI)	
ABI 7000	A
ABI 7300	A
ABI 7500 Standard 96-well Block	A
ABI 7500 FAST 96-well Block	C
ABI 7900HT Standard 96-well Block	A
ABI 7500HT FAST 96-well Block	C
ABI 7500HT 384-well Block	E
ABI 5700 (Perkin Elmer)	A
ABI 7700 (Perkin Elmer)	A
ABI StepOnePlus	C
Bio-Rad	
iCycler	A
iQ5	A
MyiQ	A
CFX96	D
Chromo4 (MJ Research)	A
Opticon (2) (MJ Research)	D
Stratagene	
Mx3005p	A
Mx3000p	A
Mx4000	D
Roche	
LightCycler 480 96-well Block	F
LightCycler 480 384-well Block	G
Eppendorf	
Mastercycler ep realplex	A

Master Mix	Catalog #
RT² SYBR[®] Green w/ ROX[®]	
2 Arrays (96-well)	PA-012
12 Arrays (96-well)	PA-012-12
24 Arrays (96-well)	PA-012-24
4 Arrays (384-well)	PA-012-8
25 mL	PA-112
RT² SYBR[®] Green w/ Fluorescein	
2 Arrays (96-well)	PA-011
12 Array (96-well)	PA-011-12
24 Arrays (96-well)	PA-011-24
4 Arrays (384-well)	PA-011-8
25 mL	PA-111
RT² SYBR[®] Green Only	
2 Arrays (96-well)	PA-010
12 Arrays (96-well)	PA-010-12
24 Arrays (96-well)	PA-010-24
4 Arrays (384-well)	PA-010-8
25 mL	PA-110

PCR Accessories	Catalog #
First Strand Kit	C-03
RNA Isolation Kit	PA-001
Nano PreAMP Kit	C-06
Nano PreAMP Primer Mixes	Various
FFPE RNA Extraction Kit	PA-023
FFPE PreAMP Kit	C-07
FFPE PreAMP Primer Mixes	Various
Primer Assays	Various

RT ² miRNA PCR ARRAYS	Human	Mouse	Size
Whole Genome (96-well)	MAH-100	MAM-100	8 Plates
Whole Genome (384-well)	MAH-3100	MAM-3100	2, 12, 24 Plates
Cancer (96-well)	MAH-102	MAM-102	2, 12, 24 Plates
Cancer (384-well)	MAH-102-E,G	MAM-102-E,G	4 Plates
Cell Differentiation/Development (96-well)	MAH-103	MAM-103	2, 12, 24 Plates
Cell Differentiation/Development (384-well)	MAH-103-E,G	MAM-103-E,G	4 Plates
miFinder™ (96-well)	MAH-001	MAM-001	2, 12, 24 Plates
miFinder (384-well)	MAH-001-E,G	MAM-001-E,G	4 Plates
RT ² qPCR-Grade miRNA Isolation Kit	MA-01	MA-01	12 Samples
RT ² miRNA First Strand Kit	MA-03	MA-03	12 Samples
RT ² miRNA qPCR Assays	Various	Various	200 Reactions

Methyl-Profiler™ DNA Methylation PCR ARRAYS	Signature Panel*	Complete Panel*
Human Breast Cancer	MeAH-011	MeAH-8010
Human Gastric Cancer	MeAH-021	MeAH-8020
Human Liver Cancer	MeAH-031	MeAH-8030
Human Lung Cancer	MeAH-041	MeAH-8040
Human Prostate Cancer	MeAH-051	MeAH-8050
Human Colon Cancer	MeAH-061	MeAH-8060
Human Stem Cell Transcription Factor	MeAH-511	
Human Inflammatory Response	MeAH-521	
Human T Cell Activation	MeAH-531	
Human Cytokine Production	MeAH-541	
Methyl-Profiler DNA Methylation Enzyme Kit (12 samples)	MeA-03	MeA-03
Human qPCR Primers for Promoter CpG Islands (200 reactions)	Inquire	Inquire
Custom Methyl-Profiler PCR Array	Inquire	Inquire

* Methyl-Profiler Signature Panels profile the promoter methylation status of 24 genes and Complete Panels profile 96 genes.

ChampionChIP™ PCR ARRAYS	Size	Catalog #
ChampionChIP PCR Arrays		
Human Stem Cell Transcription Factors	4, 12, 24 Plates	GH-501
Human Oncogene & Tumor Suppressor Genes	4, 12, 24 Plates	GH-502
Human T Helper Cell Differentiation	4, 12, 24 Plates	GH-503
ChampionChIP One-Day Kit	12 samples	GA-101
ChampionChIP Antibody Kits (antibody and qPCR controls)		
Human RNA Polymerase II	12 Samples	GA-111
Human p53	12 Samples	GA-112
Human Histone H3Ac	12 Samples	GAH-7201
Human Histone H4Ac	12 Samples	GAH-5207, GAH-3203, GAH-8208
Human Histone H3K4me1, H3K4me2, H3K4me3	12 Samples	GAH-3203
Human Histone H3K9me1, H3K9me2, H3K9me3	12 Samples	GAH-5210, GAH-8211, GAH-6204
Human Histone H3K27me3	12 Samples	GAH-9205
ChampionChIP qPCR Primers - Whole Genome	200 Reactions	Various
Custom ChampionChIP qPCR Arrays	Inquire	Inquire

FREE Data Analysis Software	Web Page
Web-based Data Analysis Software	
RT ² Profiler PCR Arrays	www.SABiosciences.com/pcrarraydataanalysis.php
RT ² miRNA PCR Arrays	www.SABiosciences.com/mirnaArrayDataAnalysis.php
EXCEL-based Data Analysis Software	
Methyl-Profiler PCR Arrays	www.SABiosciences.com/dna_methylation_data_analysis.php
ChampionChIP qPCR System	www.SABiosciences.com/chippqcr.php

Service Core for Gene Expression Analysis	Web Page
Quantitative PCR Services	
RT ² Profiler PCR Array - Focus on Pathways or Diseases	www.SABiosciences.com/PCRArrayService.php
RT ² Real-Time PCR - Analyze Individual Genes and Verify Array Data	www.SABiosciences.com/RT2PCRSERVICE.php
RNA Isolation Services	
Standard - Preparation from Cells, Tissue, Blood and Other Samples	www.SABiosciences.com/RNAIsolationService.php
Specialized - Preparation from FFPE, LCM, FNAB and Other Samples	www.SABiosciences.com/RNAIsolationService.php
Complete Service Solutions for FFPE, LCM, FNAB and Other Small Samples	
Formalin-Fixed Paraffin-Embedded	www.SABiosciences.com/FFPEService.php
Laser Capture Microdissection	www.SABiosciences.com/LCMFNABService.php
Fine Needle Aspiration Biopsies and Other Small Samples	www.SABiosciences.com/LCMFNABService.php
Illumina Gene Expression Analysis Services	
Genotyping and Copy Number Variation (CNV)	www.SABiosciences.com/illuminagenotyping.php
Whole- Genome Expression BeadChips	www.SABiosciences.com/illuminageneexp.php



ORDERING AND TECH SUPPORT

TEL: **1.888.503.3187**
1.301.682.9200

FAX: **1.888.465.9859**
1.301.682.7300

EMAIL: **Order@SABiosciences.com**
Support@SABiosciences.com

WEB: **www.SABiosciences.com**

SABiosciences Corporation
6951 Executive Way, Frederick, MD 21703 USA

THE PCR ARRAY EXPERTS™

©2009 SABiosciences Corporation All Rights Reserved.

RT² Profiler PCR Array, RT² Real-Time, RT² Nano PreAMP, RT² FFPE PreAMP, XpressRef, PathwayFinder, miFinder, Methyl-Profiler, ChampionCHIP, Focus on Your Genes, Focus on Your Pathway, The PCR Array Experts, SABiosciences, and the SA logo are trademarks of SABiosciences Corporation. ABI, ROX, and StepOnePlus are registered trademarks of Applied Biosystems. Opticon 2, Chromo4, iQ5, iCycler, CFX96, CFX384, and MyiQ are registered trademarks of BioRad Laboratories, Inc. LabChip is a registered trademark of Caliper Life Sciences. LightCycler is a registered trademark of Roche Applied Sciences. SmartCycler is a registered trademark of Cepheid. SYBR is a registered trademark of Molecular Probes. TRIzol is a registered trademark of Invitrogen. Mastercycler is a registered trademark of Eppendorf. Mx3000P, Mx3005P, and Mx4000 are registered trademarks of Stratagene. RNeasy is a registered trademark of Qiagen N.V. Actos is a registered trademark of Takeda Pharmaceutical Company, Ltd. Avandia is a registered trademark of GlaxoSmithKline. Rezulin is a registered trademark of Parke-Davis Pharmaceuticals, Ltd. BioMark is a trademark of Fluidigm Corporation.